

TritonGrip TE™ Nitrile Gloves

TritonGrip TE™ black, low-profile gloves feature exceptional tactile sensitivity. These textured nitrile gloves are suitable for all environments where barrier protection and user comfort is important. These gloves are 5 mil thick to reduce hand fatigue and enhance comfort. They are powder free, non-sterile and single use. No latex was used in the manufacturing of this product.



Specifications:

Material: Nitrile

Color: Black

Length: 9.5 in

Palm thickness: 3 mil

Finger thickness: 5 mil

Tensile strength:

- Before aging: 48 MPa

- After aging: 51 MPa

Elongation:

- Before aging: 649%

- After aging: 576%

Ordering Information:

Item #	Size	Quantity
1012-47133	Small	100/bx
1012-47134	Medium	100/bx
1012-47135	Large	100/bx
1012-47136	X-Large	100/bx

TritonGrip TE™ is manufactured in compliance with multiple international standards, including the following:

Designation	
ASTM D6319	Standard Specification for Nitrile Examination Gloves for Medical Application
ASTM D5151	Standard Test Method for Detection of Holes in Medical Gloves
ASTM F1671	Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens

Tested for use with chemotherapy drugs using ASTM

Test Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time in Minutes:
Bleomycin Sulfate (15 mg/mL)	>240
Busulfan (6 mg/mL)	>240
Carboplatin (Paraplatin) (10 mg/mL)	>240
Carmustine (3.3 mg/ml)	Not for use
Cisplatin (1.0 mg/ml)	>240
Cyclophosphamide (20.0 mg/ml)	>240
Cytarabine (100 mg/mL)	>240
Dacarbazine (10.0 mg/ml)	>240
Daunorubicin (5 mg/mL)	>240
Docetaxel (10 mg/mL)	>240
Doxorubicin Hydrochloride (2.0 mg/ml)	>240
Epirubicin (Ellence) (2 mg/mL)	>240
Etoposide (20.0 mg/ml)	>240
Fludarabine (25 mg/mL)	>240
Fluorouracil (50.0 mg/ml)	>240

Test Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time in Minutes:
Gemcitabine (38 mg/mL)	>240
Idarubicin (1 mg/mL)	>240
Ifosfamide (50.0 mg/ml)	>240
Irinotecan (20 mg/mL)	>240
Mechlorethamine (1 mg/mL)	>240
Melphalan (5 mg/mL)	>240
Methotrexate (25.0 mg/ml)	>240
Mitomycin C (0.5 mg/ml)	>240
Mitoxantrone (2.0 mg/ml)	>240
Paclitaxel (6.0 mg/ml)	>240
Rituximab (10 mg/mL)	>240
Thiotepa (10.0 mg/ml)	Not for use
Trisenox (1 mg/mL)	>240
Vincristine Sulfate (1.0 mg/ml)	>240
Fentanyl tested as follows:	
Fentanyl Citrate (100 mcg/2ml)	>240