# Care and Service of Reusable *GreenLine*® Fiber Optic Laryngoscope Blades

Not Made  
with Natural  
Rubber Latex  
For Professional  
Use Only

**CLEAN BLADE IMMEDIATELY AFTER USE TO PREVENT RESIDUES FROM DRYING**

**Intended Use:** Used for endotracheal tube intubation.

**Care and Service Document:** These instructions apply to SunMed's reusable *GreenLine*® fiber optic laryngoscope blades. SunMed has validated these directions as being capable of preparing its reusable *GreenLine*® blades for re-use. Fiber optic laryngoscope blades must be processed prior to first use, and between each use, using any of the methods presented in this care and service document:

- CLEANING + (A) COLD SOLUTION HIGH-LEVEL DISINFECTION
- CLEANING + (B) GRAVITY AUTOCLAVE STERILIZATION
- CLEANING + (C) PRE-VACUUM STEAM STERILIZATION
- CLEANING + (D) STERIS AMSCO V-PRO STERILIZATION
- CLEANING + (E) STERRAD STERILIZATION

## WARNINGS and CAUTIONS


- ⚠ DO NOT USE THE FOLLOWING: Ultrasonic clean, dry heat sterilization or Potassium Hydroxide liquids.
- ⚠ CHECK laryngoscope system functionality, and ensure blade is fully engaged on handle prior to use.
- ⚠ DISCARD any component that shows evidence of damage or deterioration.

**NOTE:** The laryngoscope blade **MUST** be removed from the handle prior to cleaning, disinfection, and/or sterilization. To clean, disinfect or sterilize the handle, consult the *Care and Service Instructions for Reusable GreenLine Fiber Optic Laryngoscope Handles with Removable Lamp Cartridge (Rev. 7)*.

## CLEANING THE BLADE PRIOR TO DISINFECTION / STERILIZATION

1. Remove laryngoscope blade from the handle head (Fig.1).
2. Rinse the blade under cool running tap water, until all visible soil is removed.
4. Immerse the blade in a presoak enzymatic cleaner solution, prepared according to manufacturer's instructions, for a minimum of two minutes.
5. After the presoak, scrub the blade using a soft bristled brush, and rinse with lukewarm running tap water for a minimum of one minute to remove any residues and visible soil.
6. Immerse the blade in an enzymatic detergent prepared according to manufacturer's instructions.
7. Scrub the blade thoroughly by using a soft bristle brush, and thoroughly rinse with cool running tap water to remove all residual detergent.
8. Dry the blade with lint-free clean cloth or filtered pressurized air.

Follow with **HIGH-LEVEL DISINFECTION** or **STERILIZATION** of blade.

## (A) COLD SOLUTION HIGH-LEVEL DISINFECTION OF THE BLADE

1. Prepare Cidex® OPA or 2.4% Glutaraldehyde solution according to manufacturer's instructions.
2. Immerse the blade in the disinfectant.
3. Thoroughly rinse by immersing the blade in purified water.
4. Dry with a lint-free clean cloth or filtered pressurized air.

Follow with testing prior to use.

## STERILIZATION OF THE BLADE

Use one of the following four sterilization methods to sterilize SunMed's reusable *GreenLine*® fiber optic laryngoscope blades.

**Steam Sterilization** Insert blade in an appropriate autoclave pouch.

|             | (B) GRAVITY DISPLACEMENT STEAM | (C) PRE-VACUUM STEAM |
|-------------|--------------------------------|----------------------|
| Temperature | 121° C (250° F)                | 132° C (270° F)      |
| Cycle Time  | 30 Minutes                     | 4 Minutes            |
| Dry Time    | 15 Minutes                     | 20 Minutes           |

## (D) STERIS® AMSCO® V-PRO® STERILIZATION SYS. (E) STERRAD® STERILIZATION SYS.

- AmSCO® V-PRO® 1 Low Temperature System
- AmSCO® V-PRO® 1 Plus Low Temperature System
- AmSCO® V-PRO® maX Temperature System
- STERRAD® 100NX® System (Std and EXPRESS cycle)
- STERRAD® NX® System (Std Cycle)
- STERRAD® 100S and 200 System (Short Cycle Outside US)
- STERRAD® 50 System

Fiber Optic blades that are compatible in STERRAD® Systems must be manufactured in or after May 2012 or have a lot number "Me" or greater. Upper case "M" corresponds to the year 2012 and lower case "e" corresponds to the manufacturing month, May. Lot number will be visible on every blade (Fig.1). Lot numbering key:

| First Letter = YEAR | Second Letter = MONTH |            |               |
|---------------------|-----------------------|------------|---------------|
| M = 2012            | a = January           | e = May    | i = September |
| N = 2013            | b = February          | f = June   | j = October   |
| O = 2014            | c = March             | g = July   | k = November  |
| P = 2015            | d = April             | h = August | l = December  |

Follow with testing prior to use.

## PRIOR TO USE

Laryngoscope systems must be tested after each reprocessing cycle and prior to use. To test:

1. Fully engage blade on compatible laryngoscope handle. If the light guide fails to light, check the batteries and lamp for correct placement. Replace with new batteries or lamp if needed.
2. If light guide fails to light following additional blade engagements, do not use laryngoscope.

Manufactured for SunMed  
2710 Northridge Dr. NW, Suite A  
Grand Rapids, MI 49544 USA  
www.Sun-Med.com

SunMed is a registered trademark in the USA  
*GreenLine* is a registered trademark of SunMed  
Made in Pakistan

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GreenLine is a registered trademark of SunMed