REUSABLE COMFORTSEAL® SILICONE MASK ASSEMBLY

Intended for use with the **AeroEclipse®**
Breath Actuated Nebulizer (BAN)

Manufactured for/Fabricado para: Monaghan Medical Corp.
5 Latour Ave., Suite 1600, Plattsburgh, NY 12901 USA
Tel. 800-833-9653 • Fax. 518-561-5088
customerservice@monaghanmed.com
www.monaghanmed.com


RX ONLY

Manufactured by/Fabricado por: Trudell Medical International
725 Baransway Drive, London, Ontario, Canada, N5V 5G4

INSTRUCTIONS FOR USE

Prior to each use, inspect the mask to make sure it is not damaged or missing parts. If you notice any damaged or missing parts, replace it immediately.

1. Remove the mask assembly from its packaging and verify the elbow adapter is firmly secured to the mask (Figure 1).
2. Ensure that the exhalation valve located at the top of the mask is intact (Figure 1).
3. Rotate the elbow adapter so it is horizontal when the mask is held in an upright position (Figure 2).
4. Remove the mouthpiece from your **AeroEclipse®** BAN and attach the elbow adapter to the mouthpiece opening as shown in Figure 3. Retain mouthpiece for future use.
5. Ensure the mask maintains a good face seal prior to operating the **AeroEclipse®** BAN.
6. Follow the operating instructions supplied with the **AeroEclipse®** BAN.
7. If you are unable to trigger the **AeroEclipse®** BAN, manually activate by depressing and releasing the green button during inspiration and releasing during expiration, or use in continuous mode.

CLEANING AND DISINFECTING INSTRUCTIONS

Healthcare providers should always use clinical judgment to select the appropriate Cleaning (1) and Disinfection (2) method and frequency, based on the patients' disease state.

1. CLEANING

- a. Remove elbow adapter from mask and wash the two parts in a solution of lukewarm clean water and liquid dish detergent, allowing the parts to soak for 15 minutes. Agitate gently.
- b. Rinse parts in lukewarm clean water.
- c. To dry, shake the parts to remove any excess water. Allow parts to air dry thoroughly before reassembly.
- d. To reassemble, insert the Elbow Adapter into the opening on the **ComfortSeal®** Mask so it is horizontal when the mask is in an upright position.

2. DISINFECTION

- Follow approved institutional procedures and processes.
- To reduce the risk of infection, clean and disinfect the mask assembly between treatments.

LOW LEVEL DISINFECTION (SINGLE PATIENT USE)

Prior to disinfection, follow Cleaning steps 1 (a) and (b) above. Do not dry or reassemble the mask.

BLEACH: Soak the two parts in 1:50 bleach to water solution (approximately 1 tablespoon bleach to 3½ cups water) for 3 minutes. Rinse well with lukewarm clean water. You may substitute a respiratory equipment disinfectant (such as Control III[†]) for the bleach solution by following the manufacturer's instructions for use. Allow to air dry thoroughly before reassembly.


ISOPROPYL ALCOHOL: Wipe the outside surfaces with a tissue dipped in 70% Isopropyl Alcohol then submerge the two parts in a soak of 70% Isopropyl Alcohol for 1 minute before shaking to remove any excess alcohol. Allow to air dry thoroughly before reassembly.

HIGH LEVEL DISINFECTION INSTRUCTIONS (FOR MULTI-PATIENT REUSE WITHIN HOSPITAL)

If the mask will be used for multiple patients, it must be cleaned and sterilized between patients.

BEFORE REPROCESSING

- Visually inspect the mask for any damaged or missing parts. If any defects are noted, dispose of the mask.
- If no defects are noted, remove the elbow adapter from the mask.

CLEANING AND DRYING

1. Rinse the mask parts under lukewarm clean running water for 30 seconds.
2. Submerge all mask parts in an enzymatic cleaning solution (such as Enzol[†]), prepared according to manufacturer's instructions. Allow masks to soak for 1 minute, then brush gently for 1 minute using a nylon bristled brush to ensure all visible soil is removed.
3. Rinse thoroughly for at least 30 seconds under clean, warm running water to remove any enzymatic cleaner residue.
4. A final rinse should be performed by submerging the device in de-ionized or sterile water for one (1) minute.
5. Shake well and allow to air dry thoroughly before proceeding to High Level Disinfection.

DANGER: Failure to properly clean the mask may result in inadequate disinfection. If mask is not visually clean, repeat the cleaning and drying steps above or dispose of the mask.

HIGH LEVEL DISINFECTION

A. LOW TEMPERATURE METHOD

- Prepare the disassembled mask prior to sterilization as described in the manufacturer's instructions. The following have been validated for high level disinfection of the **ComfortSeal®** Masks:
 - Steris[†] Hydrogen Peroxide Sterilizer: V-PRO 1, V-PRO 1 Plus, V-PRO maX, or V-PRO 60 Sterilizers
 - Reattach the elbow adapter to the mask prior to next use.

B. REVITAL-OX RESET

- Submerge the disassembled mask for 8 minutes, then change to a clean, static rinse for 1 minute.
- Allow all parts to completely dry before reassembling and reattach the elbow adapter to the mask prior to next use.

C. AUTOCLAVE (STEAM STERILIZATION)

- After cleaning and disinfection, pack the dry disassembled mask in a sterilization pouch (recommended pouch size approximately: 7"x12").

- Sterilize the disassembled mask according to the autoclave manufacturer's instructions:
 - 270°F/ 132°C (to be maintained for 4 minutes)
 - Ensure the mask components do not come in contact with other materials, metal trays or shelving in the autoclave, and that nothing is placed on top of the device during reprocessing.
 - Allow all parts to completely dry before reassembling and reattach the elbow adapter to the mask prior to next use.

PASTEURIZATION

- Secure components in mesh bag and then place in metal baskets. Each disassembled device should be placed in a separate bag.
- Pasteurize according to manufacturer's instructions. The effectiveness of this procedure was validated using an Olympic Medical[†] Pasteurmatic Compact device. Pasteur SOL II[†] enzymatic detergent was used in the 30 minute wash phase (38°C). The pasteurization phase duration was 30 minutes at a temperature of 70-72°C.
- Following pasteurization, remove components, shake off excess water and place into a sterile dryer.

POST-REPROCESSING INSPECTION

- Visually inspect the mask parts after each disinfection cycle.
- If any damaged, misshapen or seriously discolored parts are noted, the mask should be discarded.

The **ComfortSeal®** Mask may be reprocessed up to 30 times or for 6 months from date of first use (whichever is first), but then it must be replaced.

STORAGE

Keep the disinfected mask in a dry, dust-free place protected from contamination.

CAUTIONS

- Clean and disinfect/sterilize before every change in patient.
- Use of cleaning and disinfection/sterilization procedures other than those listed above have not been validated and may cause damage to the device.
- Do not disassemble the product beyond what is recommended in the cleaning instructions or damage may result.

NOTES

- This product has not been made or manufactured with BPA (Bisphenol A) or Natural Rubber Latex.
- This product contains small parts when disassembled. Do not leave it unattended with children.

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