MercuryMedical
Superior Solutions
11300 - 49th Street North
Clearwater, Florida 33762-4807 USA
www.mercurymed.com

(en)

Clear Anatomical Face Mask

Single-Patient-Use

INDICATIONS FOR USE:

The intended use of the Clear Anatomical Face Mask is to direct gas to the upper airway.

WARNINGS:

1. This mask is non-conductive and should not be used in the presence of flammable anesthetics.
2. Do not clean, soak, rinse or sterilize. To do so may result in harmful residuals or have adverse side effects on this product.
3. Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended.
4. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
5. If the patient experiences skin irritation, consult physician.
6. This Mask is not suitable for providing mechanical life support ventilation.
7. In the USA, federal law restricts this Mask to sale by, or on the order of a physician. (Rx ONLY)
8. This device is assembled with PVC components which include low levels of a phthalate as indicated. Consideration should be used when treating children or pregnant or nursing women.

CAUTIONS:

1. During use, the patient's condition must be continuously monitored.

DIRECTIONS FOR USE:

1. Connect to breathing circuit or resuscitation device with appropriate fitting.
2. Inspect all fittings to ensure secure connections.
3. If the mask contains an inflation valve, inflate or deflate the air cushion using a needleless syringe until proper patient fit is accomplished.
4. Secure by hand or use a mask strap to hold in place.
5. Dispose of according to facility protocol.

Specifications:

The MERCURY MEDICAL® Clear Anatomical Face Mask was designed to meet ISO 5356-1 requirements.

Port Size:	**Adult Large** (136mm)	22 mm I.D.
	Adult Small (127mm)	22 mm I.D.
	Child (101mm)	22 mm I.D.
	Pediatric (89mm)	22 mm I.D.
	Infant (72mm)	15 mm O.D.
	Neonatal (61mm)	15 mm O.D.

MR CONDITIONAL:

Static magnetic field of 3.0 Tesla or less
Spatial Gradient of 3000 gauss/cm or less

EC REP **EMERGO EUROPE**
Westervoortsedijk 6C
6827 AT Arnhem
The Netherlands

Medical Device · Single Patient – Multiple Use · Not Made with Natural Rubber Latex · Non Sterile · DINP · By Prescription Only · MR Conditional

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