Mercury Medical
Superior Solutions
11300 - 49th Street North
Clearwater, Florida 33762-4807 USA
www.mercurymed.com

(en)

FLOWSAFE II

Disposable **BiLevel** *CPAP* System

DESCRIPTION:
The Mercury Medical Flow-Safe II®+ device is a respiratory aid intended for use with a facemask or advanced airway and gas supplying device to elevate pressure in the patient's airway. The level of airway pressure delivered to the respiratory system may be set to increase during inspiration and decrease during expiration.

CONNECTIONS:
- Standard 22mm tubing nipple or DISS Threaded Connection
- Patient connection ISO 5356-1 – 15mm taper female and 22mm taper male

INDICATIONS FOR USE:
The Mercury Medical Flow-Safe II®+ device is intended to provide CPAP or BiLevel CPAP pressures to effectively spontaneously breathing patients in the hospital, surgery center and pre-hospital environment.

WARNINGS:
- Do not allow smoking or use unit near sparking equipment, open flame, oil or other flammable chemicals.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

CAUTIONS:
- Federal (USA) law restricts this device to sale by or on the order of a physician.
- Do not clean, soak, rinse or sterilize.
- Do not attempt to disinfect any part of this device. Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended.
- In the event of obstruction of the breathing circuit, remove the non-rebreathing valve and place on supplemental oxygen per protocol or guidelines for other accepted procedures.
- Use of Flow-Safe II+ device with non-back pressure compensated flow devices may affect input gas liter flow.
- Do not use with a vented mask **or mask with anti-asphyxia valve.** Recommend use with Mercury Medical **straight port** full face deluxe mask.
- Device not intended for long term use.
- If device becomes damaged, soiled or occluded, discard the device.
- Never wash or sterilize CPAP or BiLevel CPAP therapy. If a Flow-Safe II+ system is not immediately available, cannot be used effectively or in the event of an undesirable performance, follow oxygen therapy procedures for other accepted procedures.
- To be used only by clinicians who have been trained in BiLevel/CPAP therapy.
- Adjusting input flow rate will affect CPAP and IPAP/EPAP pressures. ALWAYS verify CPAP and IPAP/EPAP pressures with manometer.
- This device is assembled with PVC components which include low levels of a phthalate as indicated. Consideration should be used when treating children or pregnant or nursing women.
- Listed values are nominal only. Actual results may vary on device and accuracy of the flowmeter.
- This device is not intended to be used as a transport ventilator, and is not intended for life support.
- Patient SaO$_2$ should be monitored using a pulse oximeter.
- Patient EtCO$_2$ should be monitored using capnography.
- Use in CPAP mode if used in conjunction with a nebulizer.
- Recommend therapy should be started in the default CPAP mode.

CONTRAINDICATIONS:
- Respiratory arrest
- Unconscious
- Cardiogenic Shock
- Pneumothorax
- Facial Anomalies
- Facial Trauma
- Airway Obstruction

DIRECTIONS FOR PROVIDING CPAP:
Ensure device is set in CPAP mode-Green switch is set to CPAP (See Figure 1).
1. Connect O2 tubing spin nut or nipple to gas source. Spin nut connects directly to flow meter threading while nipple connects (by pushing through the spin nut) to either a Christmas Tree or O2 regulator hose barb. **(See Figure 1a.)**
2. Secure the face mask snugly to patient's face using head harness.
3. Slowly increase gas flow to reach approximately 8 cm H_2O IPAP. This is the minimum IPAP pressure required for device to function properly in the BiLevel mode. Check mask fit to patient and secure connections for leaks.
 a. Adjust the flowmeter until desired IPAP pressure is obtained. Flow of 17 LPM is required to reach maximum IPAP pressure of 12-13 cm H_2O.
 b. To decrease IPAP pressure rotate IPAP knob counterclockwise. To increase IPAP pressure rotate IPAP knob clockwise. **(See Figure 3).**
4. Effective mask seal is required for device to shift into EPAP mode from IPAP mode.
5. Dispose of according to facility protocol.

DIRECTIONS FOR PROVIDING BILEVEL:
Ensure device is set in BiLevel mode by rotating the Green switch to the BiLevel setting (See Figure 2).
1. Connect O2 tubing spin nut or nipple to gas source. Spin nut connects directly to flow meter threading while nipple connects (by pushing through the spin nut) to either a Christmas Tree or O2 regulator hose barb. **(See Figure 1a.)**
2. Secure the face mask snugly to patient's face using head harness.
3. Slowly increase gas flow to reach approximately 8 cm H_2O IPAP. This is the minimum IPAP pressure required for device to function properly in the BiLevel mode. Check mask fit to patient and secure connections for leaks.
 a. Adjust the flowmeter until desired IPAP pressure is obtained. Flow of 17 LPM is required to reach maximum IPAP pressure of 12-13 cm H_2O.
 b. To decrease IPAP pressure rotate IPAP knob counterclockwise. To increase EPAP pressure rotate EPAP knob clockwise. **(See Figure 3).**
4. Effective mask seal is required for device to shift into EPAP mode from IPAP mode.
5. Dispose of according to facility protocol.




Figure 1a.

Figure 1.

Figure 2.

Figure 3.

MEASURING PRESSURE:
- Internal pressure relief valve limits maximum CPAP pressure to 25 cm H_2O @ 25 LPM.
- Do not exceed pressure limit of manometer (25 cm H_2O).
- Manometer accuracy +/-3 cm H_2O up to 15 cm H_2O and +/- 5 cm H_2O over 15 cm H_2O.

OPERATING ENVIRONMENTAL LIMITS: -18° C to +50° C
STORAGE ENVIRONMENTAL LIMITS: -40° C to +60° C

MR CONDITIONAL:
Static magnetic field of 3.0 Tesla or less
Spatial Gradient of 3000 gauss/cm or less

UDI Unique Device Identification
Single-Patient - Multiple Use
MD Medical Device
Non Sterile
DINP
Not Made with Natural Rubber Latex
Rx Rx Prescription Only
MR Conditional

EU REP
2862
EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

*Registered U.S. Trademark of Mercury Medical Enterprises, LLC.
Protected by one or more of the following patents: US 8,522,818 B1; US 8,893,716 B2; US 9,370,635 B2; US 10,258,759 B2; US 11,383,054 B2; US 11,433,194 B2; US 12,064,551 B2; US 12,465,717 B2; Other Patent(s) Pending

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