FLOWSAFE II+®

Disposable BiLevel CPAP System

LEARN MORE

About

Flow-Safe II+® is the only available Disposable Bi-level CPAP System! The lightweight disposable feature allows for easy CPAP or Bi-level CPAP therapy set up and delivery. All devices include an integrated manometer for verifying CPAP, IPAP and EPAP pressures. The Disposable Flow-Safe II+ System is ideal for situations where backup Bi-level CPAP equipment is scarce or unavailable.

Advantages

- First Disposable Bi-level Device with CPAP or Bi-level Mode switch allows easy transition from CPAP to Bi-level Mode.
- No Extra Parts – Flow Determines CPAP or Bi-level IPAP Pressure, Easy EPAP Dial adjusts expiratory pressure.
- Verify CPAP & Bi-level Pressures with Integrated Manometer. (Available on all models)
- No Capital Equipment. (Ideal for Acute Care and Pre-hospital Settings)
- Easy to transport.



Flow-Safe II+® a solution for non-invasive mechanical ventilation equipment shortages

Disposable, Latex-Free

Flow-Safe II+® Bi-Level CPAP with Full Face Mask #10-57400 (5/Box)

EPAP Dial - When in Bi-level mode, the EPAP dial allows clinicians to set EPAP pressure

Can add a standard nebulizer (Recommended in CPAP mode only)

Built-In Manometer

Built-In Pressure Relief Valve

Green knob selects CPAP or Bi-level pressure mode

Clinical Use

- When connecting the Flow-Safe II+® Bi-level CPAP System to a flowmeter/regulator start with the green knob in the CPAP position by adjusting flow rates up and down. Various CPAP pressures can be achieved and verified with the manometer.

- For Bi-level therapy, simply switch the green knob to the Bi-level position. The CPAP pressure now becomes the Inspiratory Positive Airway Pressure (IPAP). Adjust the blue EPAP knob to obtain desired Expiratory Positive Airway Pressure.

07/2024 ML2022-443 Rev. A

Mercury Medical®
Superior Solutions

Over 55 Years of Medical Device Sales and Development

Dedicated to delivering clinically differentiated critical care technology that saves lives throughout the world.

(800) 835-6633

www.mercurymed.com