MercuryMedical Superior Technology


NIFometer®

Single-Patient-Use

Negative Inspiratory Force Meter

(en)

11300 - 49th Street North
Clearwater, Florida 33762-4807 USA
www.mercurymed.com

INDICATIONS FOR USE:

The Mercury Medical® NIFometer is used to indicate negative inspiratory force. Federal law (USA) restricts this device to sale by or on the order of a physician.

DIRECTIONS FOR USE:

1. Set NIFometer gauge to zero.
2. Connect the NIFometer to the endotracheal tube or elbow if using kit with mouthpiece.
3. To record the maximum negative inspiratory force, press down on the occlusion knob, then instruct patient to inhale. Air flow will be willfully occluded by clinician to record negative inspiratory force (see diagram).
4. Release knob immediately after measurement.
5. Memory pointer identifies and remains at the maximum inspiratory force to -40cm H_2O or to -60cm H_2O.
6. To reset memory pointer, rotate memory pointer knob counterclockwise.
7. When measurement is complete, remove NIFometer from endotracheal tube and re-connect breathing circuit per manufacturer's instructions or facility protocol.
8. Dispose of per facility protocol.

WARNINGS:

- To create and measure negative inspiratory force, air flow must be willfully occluded by clinician. Once measurement is taken, release occlusion knob to allow air flow.
- Memory pointer must be set at zero before use.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
- Do not clean, rinse, or sterilize.

CAUTIONS:

- Federal law (USA) restricts this device to sale by or on the order of a physician.
- During use, the patient's condition must be continuously monitored.
- Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended.

SPECIFICATIONS:

Range:	To -40cm H_2O or to -60cm H_2O
Variance:	± 2cm H_2O -20 to -30cm H_2O
Fitting:	15mm I.D./22mm O.D.
Weight:	30 grams nominaal
Packaged:	20/Box
	10 Kits/Box

Available as "unit only" for use with endotracheal tubes or in a kit with mouthpiece, elbow and nose clip.

Patent # US 6,854,334 B2
#83-900-0232

MercuryMedical Superior Technology

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(nl)

11300 - 49th Street North
Clearwater, FL 33762-4807 VS
www.mercurymed.com

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