Mercury Medical
Superior Solutions
11300 - 49th Street North
Clearwater, Florida 33762-4807 USA
www.mercurymed.com

PEEP Adapter

Single-Patient-Use

(en)

INDICATIONS FOR USE:

The Mercury Medical PEEP Adapter is designed to make PEEP (Positive End Expiratory Pressure) available for use on resuscitators. The adapter attaches first to the patient valve of the resuscitator and then accepts a PEEP valve having either a 22 mm male or 30 mm female connector.

DIRECTIONS FOR USE:

1. If the resuscitator has a mask attached, remove the mask.
2. Holding the patient valve in one hand, gently press and twist the adapter onto the patient valve with the other hand until the outer rubber ring of the adapter is fully seated onto the patient valve.
3. Rotate the swivel connector to assure that the center rubber bushing of the adapter is properly seated.
4. Attach PEEP valve using a twisting and pressing motion.
5. Dispose of according to facility protocol.

CAUTIONS:

- Follow PEEP valve manufacturer's instructions for PEEP application and contraindications when resuscitating a patient or administering CPAP therapy.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
- Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended.

SPECIFICATIONS:

Port Size: 22 mm Male and 30 mm Female (per ISO-5356-1)
MR Conditional: Static magnetic field of 3.0 Tesla or less
Spatial Gradient of 3000 gauss/cm or less

EC REP EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

MD Medical Device · MR Conditional · Single Patient - multiple use · Not Made with Natural Rubber Latex · Non Sterile · Rx By Prescription Only · PHT Phthalate Free

Mercury Medical
Superior Solutions
11300 - 49th Street North
Clearwater, Florida 33762-4807 USA
www.mercurymed.com

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