MercuryMedical
Gateway Business Park
11300 - 49th Street North
Clearwater, Florida 33762-4807 USA
www.mercurymed.com

Silicone Face Mask

Reusable

(en)

INDICATIONS FOR USE:

The intended use of the Silicone Face Mask is to direct gas to the upper airway.

WARNINGS:

1. This mask is non conductive and should not be used in the presence of flammable anesthetics.
2. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
3. If the patient experiences skin irritation, consult physician.
4. This Mask is not suitable for providing mechanical life support ventilation.
5. In the USA, federal law restricts this Mask to sale by, or on the order of a physician. (Rx ONLY)

CAUTIONS:

1. During use, the patient's condition must be continuously monitored.

DIRECTIONS FOR USE:

1. Connect to breathing circuit or resuscitation device with appropriate fitting.
2. Inspect all fittings to ensure connections.
3. Secure by hand to hold in place.

CLEANING INSTRUCTIONS:

1. Disconnect the mask from equipment and rinse in warm ≤ 45 °C (113 °F) water for 30 seconds.
2. Prepare Silicone rubber cleaning solution according to manufacturer's instructions.
3. A proteolytic, enzymatic detergent (e.g., Enzol) diluted in warm potable water is suggested.
4. Immerse masks in detergent for ≥ 1.0 minute ensuring all surfaces are in contact with the solution.
5. Using a clean, soft brush, mechanically remove any visible soil from the masks.
6. Rinse masks with warm water for a minimum of 30 seconds to remove contaminants or detergent residue and inspect for any residual soil.
7. Repeat cleaning steps if necessary.
8. Dry masks with clean cloth or air dry or allow to air dry for 15-20 minutes.
9. Inspect the silicone for signs of damage or discoloration. Dispose of device if damaged or worn according to facility protocol.

Note: Failure to thoroughly clean the device may result in device malfunction or cross-contamination.

STEAM STERILIZATION/HIGH LEVEL DISINFECTION:

1. Double wrap clean silicone masks in CSR sterilisation wrap.
2. Place wrapped masks on instrument tray and place in autoclave.
3. Steam sterilize masks using a gravity cycle autoclave set to 132°C (270°F) for 15 minutes followed by a 30 minute dry cycle.

NUMBER OF AUTOCLAVE CYCLES: 75

Specifications:

The MERCURY MEDICAL® Silicone Face Mask was designed to meet ISO 5356-1 requirements.

Port Size:

Mask Size #5 (126mm)	22 mm I.D.
Mask Size #4 (113.1mm)	22 mm I.D.
Mask Size #3 (97.5mm)	22 mm I.D.
Mask Size #2 (76mm)	22 mm I.D.
Mask Size #1 (76mm)	15 mm O.D.
Mask Size #1 (60mm)	15 mm O.D.
Mask Size #00 (48mm)	15 mm O.D.
Mask Size #00 (42mm)	15 mm O.D.
Mask Size #000 (35mm)	15 mm O.D.
Mask Size (43mm)	15 mm O.D.
Mask Size (37mm)	15 mm O.D.

EC REP | CE | EMERGO EUROPE | Westervoortsedijk 60 | 6827 AT Arnhem | The Netherlands | MD Medical Device | Non-Sterile | Not Made with Natural Rubber Latex | MR MRI Safe | Phthalates Free | Rx By Prescription Only

12/2023 | #83-900-0376 | Rev. 10

MercuryMedical
Gateway Business Park
11300 - 49th Street North
Clearwater, Florida 33762-4807 USA
www.mercurymed.com

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