

Pediatric Transport

Patients come in all shapes and sizes, and our smallest patients create unique challenges for Patient Transport. The Pedi Mate® Restraint Series provides a quick and efficient way to safely secure pediatric patients from 10-80 lbs to almost any cot without having to purchase additional equipment or seats.



[SHOP ONLINE](#)

FERNO® PEDI-MATE®

The Pedi Mate Pediatric Restraint System quickly adapts any ambulance cot for the safe transport of children ranging in size from 10-40 lb (4.5-18 kg).

- Fully adjustable, five-point harness system securely holds patients, providing safe restraint for transport
- Rolls compactly for convenient storage and easy deployment
- Vinyl construction is nontoxic and easy-to-clean
- Three restraint straps easily attach to any cot



[SHOP ONLINE](#)

FERNO® PEDI-MATE® PLUS

The Pedi Mate Plus includes the same features as the Pedi Mate and accommodates larger patients ranging in size from 10-100 lb (4.54-45.36 kg).

- Fully adjustable, five-point harness system securely holds patients, providing safe restraint for transport
- Rolls compactly for convenient storage and easy deployment
- Vinyl construction is nontoxic and easy-to-clean
- Three restraint straps easily attach to any cot



[SHOP ONLINE](#)

FERNO® NEO MATE

The Neo Mate Pediatric Restraint System adapts any ambulance cot to safely transport children ranging in size from 5-14 lb (2.3-6.4 kg).

- "Halo" Pad keeps infant's head stable
- Fully-adjustable, five-point harness system fits children firmly, and provides safe restraint during transport
- Three restraint straps easily attach to any cot
- Rolls up compactly for efficient storage and deployment

PARTNERS IN EMS



800.533.0523 | www.boundtree.com

As your EMS-focused distributor, Bound Tree Medical has the products, services and experience required to meet your needs. For more information, contact your Bound Tree Medical Account Manager or call 800.533.0523 or go to www.BoundTree.com