

FULL FACE CPAP & BILEVEL MASKS

The deluxe full face mask is designed with a soft double silicone mask rim that improves patient comfort while creating a good seal for reducing leaks. The elastic head harness is easy to place with hook and loop straps that adjust for better patient comfort. It includes adjustable forehead pads for additional comfort reducing pressure on the nose while improving fit. It also includes quick disconnect clips for easy removal.



Deluxe Full Face Masks w/ Straight Connector & Mask Ports

ORDERING INFORMATION

Part #	Description	Packaged
#10-57112	Large adult deluxe full face contoured mask with straight connector, 2 mask ports, head harness	5/Box
#10-57113	Medium/Small adult deluxe full face contoured mask with straight connector, 2 mask ports, head harness	5/Box
#10-57114	Small/Child deluxe full face contoured mask with straight connector, 2 mask ports, head harness	5/Box



Deluxe Full Face Masks w/ 90° Elbow, Anti-Asphyxia & Mask Ports

ORDERING INFORMATION

Part #	Description	Packaged
#10-57109	Large adult deluxe full face mask with 90° swivel port, anti-asphyxia, 2 mask ports, head harness	5/Box
#10-57110	Medium/Small adult deluxe full face mask with 90° swivel port, anti-asphyxia, 2 mask ports, head harness	5/Box
#10-57111	Small/Child deluxe full face mask with 90° swivel port, anti-asphyxia, 2 mask ports, head harness	5/Box



FULL FACE CPAP & BILEVEL MASKS

90° Anti-Asphyxia Valve Elbow

The 90° Anti-Asphyxia Valve Elbow easily and affordably converts Mercury EMS Flow-Safe[®] Deluxe Full Face CPAP Masks to Acute Care CPAP/BiLevel Masks (non-vented).

Now save approximately \$20.00 on CPAP/BiLevel Masks in the Emergency Department. Mercury Flow-Safe Deluxe Full Face Masks that arrive from EMS agencies in the ED will now convert to 90° Anti-Asphyxia CPAP/BiLevel masks (non-vented) with the new CPAP elbow.

Intended Use:

The Anti-Asphyxia Valve Elbow is an accessory to the Flow-Safe (CPAP) family of products for use with CPAP and BiLevel systems.



ORDERING INFORMATION

Part #	Description	Packaged
#10-57121	90° Non-vented Anti-Asphyxia Valve Elbow	10/Box

