



Curaplex® Defibrillator Pads

04.22



**our brand.
your advantage.**



Leads out direct connect defibrillator pads

ITEM #	DESCRIPTION	QTY
16240	ZOLL, Adult	10 pair/case
16242	Physio Control, Adult	10 pair/case
16382	Cardiac Science, Adult	10 pair/case

Multi-function defibrillator pads

ITEM #	DESCRIPTION	QTY
16384	Physio Control, Adult	10 pair/case
16383	Physio Control, Pediatric	10 pair/case
16381	ZOLL, Adult	10 pair/case
16380	ZOLL, Pediatric	10 pair/case
16251	Philips, Pediatric	10 pair/case

Features:

- Latex-free
- 60" lead wire
- Color-coded packaging
- Polymer adhesive gels provide superior contact, uniform current distribution and more effective defibrillation and pacing
- All pads have 510K for Physio Control, Zoll, Philips and HP

800.533.0523

BOUNDTREE.COM



Curaplex® Defibrillator & Electrode Compatibility

Brand of Defibrillator	Model of Defibrillator	Curaplex® Compatible Electrode Item Numbers	Description
Physio-Control	LP-15 LP-20	16242 Adult Leads Out Electrode	Adult, multifunction, radiolucent, leads out wires, pace/defibrillation, monitoring electrodes w/ Physio connectors, latex free, 60 inch lead wire and 3 year shelf life, box (10 pair)
		16384 Adult Electrode	Adult, multifunction, radiolucent, pace/defibrillation, monitoring electrodes w/ Physio connectors, latex free, Peel Pouch, 60 inch lead wire and 3 year shelf life, box (10 pair)
		16383 Pediatric Electrode	Pediatric, multifunction, radiolucent, pace/defibrillation, monitoring electrodes w/ Physio connectors, latex free, Peel Pouch, 60 inch lead wire and 3 year shelf life, box (10 pair)
	LP-1000	16242 Adult Leads Out Electrode	Adult, multifunction, radiolucent, leads out wires, pace/defibrillation, monitoring electrodes w/ Physio connectors, latex free, 60 inch lead wire and 3 year shelf life, box (10 pair)
		16384 Adult Electrode	Adult, multifunction, radiolucent, pace/defibrillation, monitoring electrodes w/ Physio connectors, latex free, Peel Pouch, 60 inch lead wire and 3 year shelf life, box (10 pair)
Zoll Medical	R-Series X-Series	16240 Adult Leads Out Electrode	Adult, multifunction, radiolucent, leads out wires, pace/defibrillation, monitoring electrodes w/ Zoll connector, latex free, 60 inch lead wire and 3 year shelf life, box (10 pair)
		16381 Adult Electrode	Adult, multifunction, radiolucent, pace/defibrillation, monitoring electrodes w/ Zoll connector, latex free, Peel Pouch, 60 inch lead wire and 3 year shelf life, box (10 pair)
		16380 Pediatric Electrode	Pediatric, multifunction, radiolucent, pace/defibrillation, monitoring electrodes w/ Zoll connector, latex free, Peel Pouch, 60 inch lead wire and 3 year shelf life, box (10 pair)
	Zoll AED PLUS	N/A	Curaplex® defibrillation electrodes are NOT compatible with this Zoll device
Cardiac Science	PowerHeart AED G3 Plus PowerHeart AED G3 Pro	16382 Adult Leads Out Electrode	Adult, multifunction, radiolucent, leads out wires w/ Cardiac Science™ connector, 48" lead wire, 30 month shelf life, latex free, box (10 pair)