



Traction Leg Splints

Relieve Nerve Tension and Swelling Pressure

The Curaplex® Traction Leg Splint is used in instances of a fracture in the femur. The traction splint is used to apply traction to the limb to reduce pain, realign the limb, and minimize vascular and neurological complication.

Key Features:

- Aluminum alloy structure
- Adjustable lengths
- Adjustable straps
- Carrying case
- Available in adult and pediatric sizes



Ordering Information:

Item #	Description	Qty
95001	Traction Leg Splint, Adult	4/bx
95002	Traction Leg Splint, Pediatric	1/bx