



sheet reference guide

| FEATURES | | 30057 | 3271-60202 | 3271-60402 | 3271-63400 | 3271-63403 | 14039MS | 3271-61202 | 3271-03490EP | 3271-62602 | 3271-63202 | 3271-63401 |
|---------------|--------------------|-------|------------|------------|------------|------------|---------|------------|--------------|------------|------------|------------|
| Size | 30" x 72" | ✓ | ✓ | ✓ | | | ✓ | ✓ | | | | |
| | 30" x 84" | | | | | ✓ | | | | | ✓ | |
| | 34" x 84" | | | | | | | | | | | ✓ |
| | 34" x 90" | | | | | | | | ✓ | | | |
| | 36" x 88" | | | | ✓ | | | | | | | |
| | 36" x 90" (XPS) | | | | | | | | | ✓ | | |
| Weight Rating | 300 lbs | ✓ | ✓ | ✓ | ✓ | ✓ | | | | | | |
| | 350 lbs | | | | | | ✓ | ✓ | | | | |
| | 400 lbs | | | | | | | | ✓ | ✓ | ✓ | ✓ |
| Material | Fluid Resistant | ✓ | ✓ | ✓ | ✓ | ✓ | | ✓ | ✓ | ✓ | ✓ | ✓ |
| | Fluid Impervious | | | | | | ✓ | | | | | |
| End Type | Elastic ends | | ✓ | | ✓ | ✓ | | ✓ | ✓ | ✓ | ✓ | ✓ |
| | Square ends | ✓ | | ✓ | | | ✓ | | | | | |
| Color | Blue | ✓ | ✓ | ✓ | ✓ | ✓ | | | ✓ | ✓ | ✓ | ✓ |
| | White | | | | | | ✓ | ✓ | | | | |

