




Ambu® SureSight® Mobile for pre-hospital settings

EASY-TO-USE INTERFACE

SureSight Mobile's intuitive design can help shorten the learning curve and reduce task complexity, which can be crucial for patient safety in unexpected and emergent cases.

- **Familiar and intuitive:** Ergonomic design and simple touch screen experience.
- **Capture images and videos at the push of a button** to document intubation.
- **Verify your intubation:** Simple playback feature allows recording of intubation to be reviewed without needing to transfer files.

THE GO-TO VIDEO LARYNGOSCOPE SOLUTION FOR EMERGENCY MEDICINE PROFESSIONALS

SureSight Mobile is a handheld video laryngoscope that offers flexibility, reliability, and a positive user experience, all of which helps you intubate with confidence in a variety of pre-hospital settings.

A VIDEO LARYNGOSCOPE FOR A VARIETY OF PRE-HOSPITAL SETTINGS

Portability, versatility, and a small footprint, combined with high-quality imaging, make SureSight Mobile an ideal option for emergent and unplanned intubations.

- **Optimize the image for different environments:** Brightness controls and tiltable screen enable you to enhance the image whether you are in the ambulance or in full sunlight.
- **Convenient carry case:** A portable storage solution with a small footprint for environments with limited storage space.
- **The comprehensive range of blades** for adult and pediatric patients is compatible with both the SureSight Connect and SureSight Mobile solutions.

ROBUST AND RELIABLE

SureSight Mobile is designed to be robust and reliable for use under challenging conditions.

- **Designed for durability** in a wide range of environments and terrains.
- **Audible click confirms correct attachment** of blade to SureSight Mobile display.
- **High-quality imaging** helps you clearly visualize the airway anatomy.
- **Continuous usage:** Exchangeable and rechargeable batteries with up to 100 minutes of battery operation time per battery that minimizes downtime risk.

Questions? Contact your Ambu representative at **800-262-8462** or visit [AmbuUSA.com/SureSight](https://www.AmbuUSA.com/SureSight)

SPECIFICATIONS

Display	
Display Type	3.5" TFT Touch LCD
Startup Time	3s
Memory	
Storage Capacity	64GB
Electrical System (SureSight Battery)	
Battery Type	Li-Ion 3.6 Vdc 2200 mAh
Battery Operation Time	100 minutes
Electrical System (SureSight Battery Charger)	
Power Requirement	5.0 Vdc, 3.0 A, 15.0 W
SureSight Data Cable Specification	
Cable Length	1.2 meters (3.93 ft)
Transport Conditions	
Temperature [°C (°F)]	-10 - 55 (14 - 131)
Relative Humidity [%]	10 - 95
Atmospheric Pressure [kPa, (psi)]	50 - 106 (7 - 15)
Storage Conditions	
Temperature [°C (°F)]	10 - 25 (50 - 77)
Relative Humidity [%]	10 - 85
Atmospheric Pressure [kPa, (psi)]	50 - 106 (7 - 15)
Operating Environment	
Temperature [°C (°F)]	0 - 40 (32 - 104)
Relative Humidity [%]	15 - 95
Atmospheric Pressure [kPa, (psi)]	62 - 106 (9 - 15)
Altitude [m]	≤ 4,000
Operating Environment (SureSight Battery Charger and Data Cable)	
Temperature [°C (°F)]	10 - 40 (50 - 104)
Relative Humidity [%]	30 - 85
Atmospheric Pressure [kPa, (psi)]	80 - 106 (12 - 15)
Altitude [m]	≤ 2,000

ORDERING SPECIFICATIONS

Item No.	Product Description	Quantity/Box
390471000	SureSight Mobile Kit - US (SureSight Mobile w. Battery, Carry Case, Extra Battery, Battery Charger and Data Cable)	1
390481000	SureSight Mobile Mini Kit - US (SureSight Mobile w. Battery, Carry Case and Extra Battery)	1
390611000	SureSight Battery - US	1
390621000	SureSight Battery Charger - US	1
390631000	SureSight Data Cable - US	1
390641000	SureSight Carry Case - US	1
390101000	Ambu SureSight Blade Miller Size 0 US	5
390111000	Ambu SureSight Blade Miller Size 1 US	5
390211000	Ambu SureSight Blade Hyperangulated Size 1 US	5
390221000	Ambu SureSight Blade Hyperangulated Size 2 US	5
390521000	Ambu SureSight Blade Hyperangulated Size 2 Channeled US	5
390231000	Ambu SureSight Blade Hyperangulated Size 3 US	10
390531000	Ambu SureSight Blade Hyperangulated Size 3 Channeled US	10
390241000	Ambu SureSight Blade Hyperangulated Size 4 US	10
390331000	Ambu SureSight Blade MAC Size 3 US	10
390341000	Ambu SureSight Blade MAC Size 4 US	10



SureSight Mobile with Hyper 3 Channeled blade eliminates the need for a stylus



SureSight Mobile features 200° touch-screen adjustment

Ambu

Ambu, Inc.
6721 Columbia Gateway Drive Suite 200
Columbia, Maryland 21046
Tel. 800 262 8462
Fax 800 262 8673
ambuUSA.com



The blades are PC bioplastic material (using bio-attribution) and have 55% second-generation bio-based feedstock and 45% traditional fossil-based feedstock (certified under mass-balance).