

# Air-Q®sp3G Intubating Laryngeal Airway with Gastric Access

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

## Device Manufacturing Specifications

Reference Number	60005, 60055, 60105, 60205, 60305, 60405, 60505
Manufacturer	AirLife
Classification - US	Class I
FDA Product Code	CAE-Airway, Oropharyngeal, Anesthesiology
Classification - EU	Class IIa
CE Mark/Notified body	CE0482/Medcert GmbH
EU Authorized Representative	Mt Promedt Consulting GmbH
GMDN Code	42424 Oropharyngeal Airway, Single Use
UMDNS Code	10059 Airways, Oropharyngeal
Made In	China
Usage	Disposable, Single Patient Use
Sterile	EO Sterilization
Patient Population	Infant, Pediatric, Adult
Packaging	Individually Packaged, 10/Case

**Description:** A curved tube used in inhalational anesthesia and resuscitation to facilitate and secure airway patency for the delivery and exchange of gases in spontaneously breathing and ventilated patients. Air-Q®sp3G has a self pressurizing cuff and side channels for gastric access.

**Intended Purpose:** To provide an unobstructed airway for gas delivery.

**Area of Use:** Hospitals, pre-hospital, surgical centers.

**Contraindications:** In patients at high risk for regurgitation and/or aspiration.

## Device Specifications

Size	Adapter Color	Ideal Body Weight	Maximum OEETT	Mouth Opening	Length of Tube	Maximum OG Tube	ET Tube Connector
0	Light Purple	<2.0kb	3.0 mm	5 mm	6.0 cm	5 Fr	15 mm
0.5	Pink	2 – 4 kg	4.0 mm	8 mm	7.0 cm	6 Fr	15 mm
1.0	Blue	4 – 7 kg	4.5 mm	11 mm	9.0 cm	8 Fr	15 mm
1.5	Green	7 – 17 kg	5.0 mm	14 mm	11.0 cm	8 Fr	15 mm
2	Orange	17 – 30 kg	5.5 mm	17 mm	14.0 cm	10 Fr	15 mm
3	Yellow	30 – 60 kg	7.0 mm	20 mm	16.0 cm	14 Fr	15 mm
4	Red	60 – 80 kg	8.0 mm	23 mm	18.0 cm	16 Fr	15 mm
5	Purple	> 80 kg	9.0 mm	25 mm	20.0 cm	18 Fr	15 mm



## Device Material

Component	Material
Airway Connector	Polycarbonate
Mask	Silicone
Airway Tube	Silicone

**Latex:** AirLife® does not utilize latex or latex-containing material to manufacture products. To the best of our knowledge latex does not contact components or finished goods during the manufacturing process.

**Phthalates:** The selected materials are Non-DEHP; DEHP was not intentionally added as part of the manufacturing process.

**Mercury-Lead:** The selected materials do not contain lead or mercury.

**Biocompatible:** Per device classification in ISO 10993. AirLife manufactures product from medical grade materials in compliance with Good Manufacturing Practices (GMP's) as listed in 21 C.F.R. (U.S. code of Federal Regulations).

Part Number	Size	UOM	GTIN	UOM	Case GTIN
2114-60005	0	Each	00814954020263	Case	10814954020260
2114-60055	0.5	Each	00814954020270	Case	10814954020277
2114-60105	1.0	Each	00814954020287	Case	10814954021083
2114-60155	1.5	Each	00814954020294	Case	10814954021090
2114-60205	2	Each	00814954020300	Case	10814954021106
2114-60305	3	Each	00814954020317	Case	10814954021113
2114-60405	4	Each	00814954020324	Case	10814954021120
2114-60505	5	Each	00814954020331	Case	10814954021137