

EtCO₂ sampling cannulas


EtCO₂/O₂ sampling cannulas collect the exhaled CO₂ while simultaneously delivering oxygen. These EtCO₂/O₂ cannulas also feature a reflective style connector, compatible with Microstream™, Capnostream™ and Oridion™ capnography monitors.

- Simultaneously deliver oxygen while obtaining CO₂ sampling during spontaneous breathing
- Soft cannula and tubing material help to ensure patient comfort
- 3-channel tubing prevents kinking
- Tubing length of 7', with Fits-All Universal O₂ connector
- In-line filter for improved moisture management

sampling line with airway adapter (T-piece)

Sampling line with Airway Adapter (T-Piece) can be used on both intubated and non-intubated patients. Available in standard use. Can be used on Adult/Pediatric patients with ET tube size > 4.5 mm.

In-line hydrophobic



In-line hydrophobic filter for improved moisture management

Orange reflective connector indicates the standard-use cannulas

Reflective style connector with ergonomic grip

Curved nares made of ultra-soft material

Fits-All O₂ connector comes standard

The tubing and ear piece are designed to help prevent kinking

Soft-Ears™ material helps ensure superior patient comfort, even during prolonged use

Item #	Description	Qty
301-5707F-SE	Adult EtCO ₂ /O ₂ Oral/Nasal Cannula, Dual Sampling	25/ca
301-5207F-SE	Pediatric EtCO ₂ /O ₂ Nasal Cannula, Dual Sampling	25/ca
301-5807F-SE	Pediatric EtCO ₂ /O ₂ Oral/Nasal Cannula	25/ca
301-5107ET	Adult/Pediatric Sampling Line w/ Airway Adapter Intubated	25/ca

Kits		
670270-KIT	Curaplex® EtCO ₂ /O ₂ Oral-Nasal Cannula with Universal Connector - 5pk Multipack	30/ca
670272-KIT	Curaplex® EtCO ₂ /O ₂ Oral/Nasal Cannula 10pk Multipack	15/ca
670271-KIT	Curaplex® EtCO ₂ Sampling Line with Airway Adapter - 5pk Multipack	20/ca
670273-KIT	Curaplex® Sampling Lines with Airway Adapter 10pk Multipack	25/ca

Microstream™, Capnostream™, and Oridion™ are trademarks of Medtronic, related to Product marked by Medtronic.