



Sheets

Our sheets are strong enough to assist in patient transfers of various weights and are resistant to moisture and fluids.

Features:

- Durable and cost-effective
- Designed to transport patients and protect equipment
- Made from spunbound polypropylene
- Fitted and flat sheets, pillow cases, and sets available
- Multiple sizes and weight ratings from 300 - 400 lbs.

Fitted Sheets

Item #	Weight Rating (lbs.)	Description
3271-62602	400	Fitted Sheet, Stryker XPS, 36 in x 90 in, Elastic Ends, Blue, 30/CS
3271-63202	400	Fitted Sheet, 30 in x 84 in, Elastic Ends, Blue 50/CS
3271-63401	400	Fitted Sheet, 34 in x 84 in, Elastic Ends, Blue, 50/CS
3271-63402	400	Fitted Sheet, 30 in x 84 in, Elastic Ends, Blue, Fluid Impervious, 50/CS
3271-03490EP	400	Fitted Sheet, Ferno® PowerFlexx, 34 in x 90 in, Elastic Ends, Blue, 50/CS
14039MS	350	Fitted Sheet, 30 in x 72 in, Square Ends, White, 50/CS
3271-61202	350	Fitted Sheet, 30 in x 72 in, Elastic Ends, White, 50/CS
3271-63400	300	Fitted Sheet, 36 in x 88 in, Elastic Ends, Blue, 50/CS
3271-63403	300	Fitted Sheet, 30 in x 84 in, Elastic Ends, Blue, 50/CS
3271-60402	300	Fitted Sheet, 30 in x 72 in, Square Ends, Blue, 50/CS
3271-60202	300	Fitted Sheet, 30 in x 72 in, Elastic Ends, Blue, 50/CS
30057	300	Fitted Sheet, 30 in x 72 in, Square Ends, Blue, 50/CS



3271-62602

Curaplex 36 in x 90 in
Stryker XPS Fitted Sheets

Pillow Cases


Item #	Description
3271-44125	Pillow Case, 30 in x 22 in, White, Poly/Tissue
3271-44150EP	Pillow Case, 30 in x 22 in, White, 100/CS

Linen Sets/Kits

Item #	Description
3271-300LSEP	3-Piece Linen Set w/ Fitted Sheet, Flat Sheet, and Pillow Case, 25/CS
8600-01339	3-Piece Linen Set w/ Stryker XPS Fitted Sheet, Flat Sheet, and Pillow Case, 30/CS
670193-KIT	3-Piece Bedding Kit w/ Fitted Sheet, Flat Sheet, and Pillow Case, 50/CS


Flat Sheets

Item #	Description
3271-44100	Flat Sheet, 40 in x 90 in, Blue, Poly/Tissue, 50/CS
3271-61002	Flat Sheet, 50 in x 84 in, Navy Blue, 50/CS



Sheet Reference Guide

FEATURES		30057	3271-60202	3271-60402	3271-63400	3271-63403	14039MS	3271-61202	3271-03490EP	3271-62602	3271-63202	3271-63401	3271-63402
Size	30 in x 72 in	✓	✓	✓			✓	✓					
	30 in x 84 in					✓					✓		✓
	34 in x 84 in											✓	
	34 in x 90 in (PowerFlexx)								✓				
	36 in x 88 in				✓								
	36 in x 90 in (XPS)									✓			
Weight Rating	300 lbs	✓	✓	✓	✓	✓							
	350 lbs						✓	✓					
	400 lbs								✓	✓	✓	✓	✓
Material	Fluid Resistant	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
	Fluid Impervious												✓
End Type	Elastic ends		✓		✓	✓		✓	✓	✓	✓	✓	✓
	Square ends	✓		✓			✓						
Color	Blue	✓	✓	✓	✓	✓			✓	✓	✓	✓	✓
	White						✓	✓					



Curaplex Fitted Sheets - Stretcher/Cot Cross Reference

Multiple weight ratings are available by size (300 / 350 / 400 lbs).

Manufacturer	Model	Curaplex SKU - Square	Sheet Size	Curaplex SKU - Elastic	Sheet Size
Stryker	XPS			3271-62602	36" x 90"
Stryker	Power-PRO XT (6500/6506)	3271-60402, 3271-61202, 30057, 14039MS	30" x 72"	3271-63202, 3271-63402, 3271-63403	30" x 84"
Stryker	Power-PRO 2 (6507)	3271-60402, 3271-61202, 30057, 14039MS	30" x 72"	3271-63202, 3271-63402, 3271-63403	30" x 84"
Stryker	Performance-PRO XT (6085/6086)	3271-60402, 3271-61202, 30057, 14039MS	30" x 72"	3271-63202, 3271-63402, 3271-63403	30" x 84"
Stryker	MX-PRO R3 (6082)	3271-60402, 3271-61202, 30057, 14039MS	30" x 72"	3271-63202, 3271-63402, 3271-63403	30" x 84"
Stryker	MX-PRO Bariatric (6083)			3271-03490EP	34" x 90"
Ferno	POWERFlexx			3271-03490EP	34" x 90"
Ferno	POWER X1 / POWER X2			3271-63202, 3271-63402, 3271-63403, 3271-03490EP	30" x 84" 34" x 90"
Ferno	93-H PROFlexx			3271-63202, 3271-63402, 3271-63403, 3271-03490EP	30" x 84" 34" x 90"
Ferno	35-X PROFlexx			3271-63202, 3271-63402, 3271-63403, 3271-03490EP	30" x 84" 34" x 90"
Ferno	35-A Mobile Transporter			3271-63202, 3271-63402, 3271-63403	30" x 84"
Ferno	iNX Integrated Transport	3271-60402, 30057, 14039MS, 3271-61202	30" x 72"	3271-63202, 3271-63402, 3271-63403	30" x 84"



boundtree.com | 800.533.0523
or contact your Bound Tree account manager

- XPS = Expandable Patient Surface (Stryker). Select the 36" x 90" XPS sheet (3271-62602) ONLY when the cot is XPS-equipped.
- 3271-03490EP is the Ferno® PowerFlexx 34" x 90" sheet; it also fits several Ferno PROFlexx / POWER models as an alternate.