

CURAPLEX® ACCUCUFF™



AccuCuff™ AccuCuff Cuff Pressure Indicator

Representing a significant advance in cuff pressure monitoring, the AccuCuff™ cuff pressure indicator is one of our most innovative intubation safety devices.

Features:

- Continuous visual display eliminates the need for constant spot-checking and reduces the risk of Post Ventilatory Pneumonia
- Negative pressure graduation ensures cuff is fully deflated, reducing the risk of trauma during removal
- Lightweight device reduces drag and potential tube displacement
- Suitable for use on a wide range of airway devices extending usage into many clinical areas
- Simple design for ease of use
- Pressure can be quickly and easily adjusted avoiding excessive cuff inflation
- Individually wrapped
- MRI compatible
- Single use
- Latex free
- 5 year shelf life allows long term storage and reduces waste
- Supplied in boxes of 50

Item #	Description	Qty
301-AC0100B-BX	Curaplex® AccuCuff™ Tracheal Cuff Pressure Indicator, Luer Lock (Blue)	50/bx



Contact your local Tri-anim Account Manager for more information or visit: tri-anim.com