

TritonGrip SE™ Nitrile Gloves



TritonGrip SE™ nitrile gloves feature exceptional tactile sensitivity. Suitable for environments where barrier protection and user comfort is important, these flexible, textured gloves are made to be comfortable all day. They are powder free, non-sterile and single use. No latex was used in the manufacturing of this product.

Ordering Information

1012-47123	Small, 100/bx
1012-47124	Medium, 100/bx
1012-47125	Large, 100/bx
1012-47126	X-Large, 90/bx

Specifications

Material: Nitrile
Color: Blue

Length: 9.5 in
Cuff thickness: 5 mil
Palm thickness: 5 mil
Finger thickness: 6 mil

Tensile strength
- Before aging: 21 MPa min
- After aging: 21 MPa min

Elongation
- Before aging: 550%
- After aging: 500%

For more information, please contact your Account Manager or call Customer Service at 800.533.0523

