

TritonGrip TE™ Nitrile Gloves



TritonGripTE™ black, low-profile gloves feature exceptional tactile sensitivity. These textured nitrile gloves are suitable for all environments where barrier protection and user comfort is important. These gloves are 5 mil thick to reduce hand fatigue and enhance comfort. They are powder free, non-sterile and single use. No latex was used in the manufacturing of this product.

Ordering Information

1012-47133	Small, 100/bx
1012-47134	Medium, 100/bx
1012-47135	Large, 100/bx
1012-47136	X-Large, 90/bx

Specifications

Material: Nitrile
Color: Black

Length: 9.5 in
Cuff thickness: 5 mil
Palm thickness: 4 mil
Finger thickness: 5 mil

Tensile strength
- Before aging: 35 MPa min
- After aging: 38 MPa min

Elongation
- Before aging: 580%
- After aging: 550%

For more information, please contact your Account Manager or call Customer Service at 800.533.0523

