

AnsellBound Tree

# Responding to the needs of **EMS PROFESSIONALS**

## **MICRO-TOUCH® Royal Blue Nitrile is transitioning to MICROFLEX® UltraSense™ US-220**

EMS professionals are constantly in need of reliable protection to provide safety to others. They face hazards that could be life threatening, which is why there should be no compromise on quality personal protective equipment (PPE). Our transition from MICRO-TOUCH® Royal Blue Nitrile to MICROFLEX® UltraSense™ US-220 comes with the same levels of protection and comfort that you can rely on.



**SAME RELIABLE  
PROTECTION,  
NEW LOOK**

➤ To learn more, visit [Ansell.com/microflex-ultrasense-us-220](https://Ansell.com/microflex-ultrasense-us-220)

**MICROFLEX®**



# Ansell

## THE SAME PROTECTION YOU CAN RELY ON

MICROFLEX® UltraSense™ is a great first choice for mid-duty level protection. They are also designed to resist sticking to tape and adhesives.

### MICROFLEX®



#### MICROFLEX® UltraSense™ US-220

Designed to provide reliable protection and tactility

Standard

0.08mm / 3.1mil

0.13mm / 5.1mil

Textured Fingers

- Designed with a thin, yet durable nitrile formulation for dependable protection
- Royal blue color provides greater contrast than standard light blue gloves
- Non-stick exterior surface prevents tapes and adhesives from sticking
- Approved for use with chemotherapy drugs

Nitrile

1.5 AQL

Yes, in accordance with ASTM D6978 and US FDA cleared

Static Testing used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes)  
ACPP (Dynamic Testing Protocol)\* used for Fentanyl testing - Results showed no permeation after 30 min. of dynamic testing protocol  
\*Ansell Cytostatic Permeation Program (ACPP) - ACPP testing was designed to mimic in-use working conditions and worst case scenarios providing glove users a perspective for safer handling of harmful or hazardous drugs

US QSR/FDA 510k

Compliant

Not Tested

XS (5.5 - 6): US220XS L (8.5 - 9): US220L  
S (6.5 - 7): US220S XL (9.5 - 10): US220XL  
M (7.5 - 8): US220M

100 gloves per dispenser / 10 dispensers per case / 1000 gloves per case

### MICRO-TOUCH®



#### MICRO-TOUCH® ROYAL BLUE NITRILE

Designed to provide reliable protection and tactility

Standard

0.08mm / 3.1mil

0.13mm / 5.1mil

Textured Fingers

- Designed with a thin, yet durable nitrile formulation
- Royal blue color provides greater contrast than standard light blue nitrile gloves
- Approved for use with chemotherapy drugs

Nitrile

1.5 AQL

Yes, in accordance with ASTM D6978 and US FDA cleared

Static Testing used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes)  
ACPP (Dynamic Testing Protocol)\* used for Fentanyl testing - Results showed no permeation after 30 min. of dynamic testing protocol  
\*Ansell Cytostatic Permeation Program (ACPP) - ACPP testing was designed to mimic in-use working conditions and worst case scenarios providing glove users a perspective for safer handling of harmful or hazardous drugs

US QSR/FDA 510k

Compliant

Not Tested

S (6.5 - 7): 1015-29070 L (8.5 - 9): 1015-29090  
M (7.5 - 8): 1015-29080 XL (9.5 - 10): 1015-29100

100 gloves per dispenser / 10 dispensers per case / 1000 gloves per case

Boundtree.com

Contact your Boundtree representative for ordering or more information.

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates.  
US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2022 Ansell Limited. All Rights Reserved.

Bound Tree