

Model S-240



Model S-242



FEATURES AND BENEFITS

The **Model S-240, S-242** is a compact, foldable transport chair designed for efficient storage in tight spaces while offering secure and reliable patient handling. It features telescoping front handles with two-position adjustments and rear handles—two on the S-240 and four on the S-242—for improved control and maneuverability. The S-240 model includes two 5-inch rear wheels, while the S-242 adds two 3-inch locking wheels and two 5-inch wheels for enhanced mobility and stability. The durable, vinyl-coated nylon cover resists stains, blood, and bodily fluids and is easy to clean. Equipped with three quick-release patient restraints and leg restraints, the chair ensures safe and secure transport.

SPECIFICATIONS

Item Number S-240	0721324
Item Number S-242	0721323
Length	29 in / 740 mm
Length (folded)	36 in / 920 mm
Width	21 in / 530 mm
Width (folded)	21 in / 530 mm
Height	36 in / 920 mm
Height (folded)	7 in / 170 mm
Weight	19 lb / 8,65 kg
Weight S-240 (Without restraints)	19 lb / 8,65 kg
Weight S-242 (Without lower rear handles)	22 lb / 10 kg
Load Limit	375 lb / 170 kg