
Discover comprehensive solutions for pleural and pericardial drainage









Contents

PLEURAL DRAINAGE CATHETERS

Pneumothorax	
Wayne Pneumothorax Catheter Set and Tray – Seldinger, Pigtail	4
Wayne Pneumothorax Catheter Set – Trocar, Pigtail	5
Cook® Emergency Pneumothorax Set	6
Pneumothorax Set and Tray	7
Catheter Aspiration Set for Simple Pneumothorax	8
Fluid drainage	
Thal-Quick Chest Tube Set	10
Thal-Quick Chest Tube Tray	11
MULTIPURPOSE CATHETERS	
Fuhrman Pleural/Pneumopericardial Drainage Set – Seldinger, Pigtail	12
Fuhrman Pleural/Pneumopericardial Drainage Tray – Seldinger, Pigtail	13
PERICARDIOCENTESIS CATHETERS	
Lock Pericardiocentesis Set and Tray – Pigtail Catheter	14
PERITONEAL LAVAGE CATHETERS	
Peritoneal Lavage Set and Tray	15
ACCESSORIES	
Simple Pneumothorax Aspiration Accessory Set	16
Plastic Multipurpose Tubing Adapter	17
Vinyl Connecting Tube	17
Thal-Quick Chest Tube Adapter	18



Wayne Pneumothorax Catheter Set and Tray – Seldinger, Pigtail


Used for the relief of simple, spontaneous, iatrogenic, and tension pneumothorax.

Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The set comes with the Cook Chest Drain Valve, which is designed to help remove air from the pleural cavity.
- The enlarged sideports are designed to facilitate the drainage of air.
- The radiopaque catheter material enhances x-ray visualization.
- Pigtail catheters have a comparable efficacy to chest tubes in patients with pneumothorax.¹
- The catheter comes with a pre-attached, three-way clear plastic stopcock.
 - 1. Kulvatunyou N, Vijayasekaran A, Hansen A, et al. Two-year experience of using pigtail catheters to treat traumatic pneumothorax: a changing trend. *J Trauma*. 2011;71(5):1104-1107.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr
Seldinger Set									
G56535	C-UTPT-1400-WAYNE-112497-IMH	14	29	19	18	5, 9	0.038	70	14
G56536	C-UTPT-1400-WAYNE-A-112497-IMH*	14	29	19	18	5	0.038	70	14
Seldinger Tra	у								
G56537	C-UTPTY-1400-WAYNE-112497-IMH	14	29	19	18	5, 9	0.038	70	14

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Support & Distribution for details. *This set does not include the Cook Chest Drain Valve.





Wayne Pneumothorax Catheter Set – Trocar, Pigtail

Used for the relief of simple, spontaneous, iatrogenic, and tension pneumothorax.


Features and benefits

- The needle obturator comes pre-assembled within the Wayne Pneumothorax Catheter to allow for direct-stick insertion.
- The set comes with the Cook Chest Drain Valve, which is designed to help remove air from the pleural cavity.
- The enlarged sideports are designed to facilitate the drainage of air.
- The radiopaque catheter material enhances x-ray visualization.
- Pigtail catheters have a comparable efficacy to chest tubes in patients with pneumothorax.*
- The catheter comes with a pre-attached, three-way clear plastic stopcock.

*Kulvatunyou N, Vijayasekaran A, Hansen A, et al. Two-year experience of using pigtail catheters to treat traumatic pneumothorax: a changing trend. *J Trauma*. 2011;71(5):1104-1107.

Trocar

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Trocar Needle gage	Needle Length cm	
G56532	C-UTPT-1020-WAYNE-IMH	10.2	29	20	19	16.5	
G56533	C-UTPT-1400-WAYNE-IMH	14	29	20	19	16.5	





Cook® Emergency Pneumothorax Set


Used for emergency relief and temporary management of suspected tension pneumothorax.

Features and benefits

• The set comes with the Cook Chest Drain Valve, which is designed to help remove air from the pleural cavity.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Needle gage	
Set with BD Ch	loraPrep Hi-Lite Orange®				
G13722	C-TPTS-8.5-6.0-FSNS	8.5	6	15	
G55524	C-TPTS-8.5-8.7-FSNS	8.5	8.7	15	
Dry Set					
G27292	C-TPTSJ-8.5-6.0-FSNS	8.5	6	15	
G55525	C-TPTSJ-8.5-8.7-FSNS	8.5	8.7	15	
Pre-Hospital Se	et*				
G29467	C-TPTS-8.5-6.0-PHS	8.5	6	15	

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Support & Distribution for details. *The set is supplied in a slip-peel pouch for transportation in emergency vehicles.



Hi-Lite Orange is a registered trademark of CareFusion 2200, Inc. Ty-Rap is a registered trademark of ASEA Brown Boveri Ltd.



Pneumothorax Set and Tray

Used for the relief of simple, spontaneous, or iatrogenic pneumothorax and tension pneumothorax.


Features and benefits

- The straight catheter is used for trocar placement.
- The set comes with the Cook Chest Drain Valve, which is designed to help remove air from the pleural cavity.
- The radiopaque catheter material enhances x-ray visualization.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	
Set							
G03465	C-TPTP-100	6.3	19	6	21	12	
G03301	C-TPT-100	9	29	20	18	20	
G03302	C-TPT-200	10.2	29	20	18	20	
Tray							
G03311	C-TPTY-100	9	29	20	18	20	

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Support & Distribution for details.

Set components





Catheter Aspiration Set for Simple Pneumothorax


Used for the relief of simple, spontaneous, or iatrogenic pneumothorax and tension pneumothorax.

Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The radiopaque catheter material enhances x-ray visualization.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Wire Guide Diameter inch	Wire Guide Length cm	Needle gage	Needle Length cm
G05296	C-CASP-901-FORD	9	25	10	0.038	50	18	7

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Support & Distribution for details.





THREE-WAY CLEAR PLASTIC STOPCOCK



Thal-Quick Chest Tube Set


Used for the percutaneous introduction of a chest tube for pleural fluid drainage.

Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The centimeter-marked dilators allow for additional control over the insertion of a chest tube into the pleural space.*
- The sideports are positioned at the distal end of the chest tube.
- The clear polyvinylchloride catheter with a radiopaque stripe allows fluid to be seen during the drainage procedure.
- The most proximal side hole is positioned on the radiopaque stripe to aid in confirming catheter positioning via x-ray.
- The double lumen set has an additional lumen available for aspiration or infusion.
- The Thal-Quick Chest Tube Adapter is sold separately (see page 18).

^{*}Double lumen dilators are not marked.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr
Single Lumen Se	et								
G06786	C-TQTS-800	8	18	3	18	5	0.038	55	9
G06787	C-TQTS-1000	10	20	3	18	5	0.038	55	9, 11
G05161	C-TQTS-1200	12	22	3	18	7	0.038	55	9, 13
G05777	C-TQTS-1400	14	41	4	18	7	0.038	65	9, 16
G04891	C-TQTS-1600	16	41	4	18	7	0.038	65	9, 13, 18
G05453	C-TQTS-1800	18	41	4	18	7	0.038	65	9, 14, 20
G05454	C-TQTS-2000	20	41	4	18	7	0.038	65	9, 14, 22
G04137	C-TQTS-2400	24	41	4	18	7	0.038	65	9, 18, 26
G04139	C-TQTS-2800	28	41	4	18	7	0.038	65	11, 20, 30
G07243	C-TQTS-3200	32	41	4	18	7	0.038	65	11, 22, 34
Double Lumen S	Set								
G36224	C-DLCTS-2400	24	40	4	18	7	0.038	80	12, 18, 26






Thal-Quick Chest Tube Tray

Used for the percutaneous introduction of a chest tube for pleural fluid drainage.

Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The dilators allow for additional control over the insertion of a chest tube into the pleural space.
- The sideports are positioned at the distal end of the chest tube.
- The clear polyvinylchloride catheter with a radiopaque stripe allows fluid to be seen during the drainage procedure.
- The most proximal side hole is positioned on the radiopaque stripe to aid in confirming catheter positioning via x-ray.
- The double lumen set has an additional lumen available for aspiration or infusion.
- The Thal-Quick Chest Tube Adapter is sold separately (see page 18).

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr
G06791	C-TQTSY-800	8	18	3	18	5	0.038	55	9
G06885	C-TQTSY-1000	10	20	3	18	5	0.038	55	11
G05464	C-TQTSY-1200	12	22	3	18	5	0.038	70	6, 10, 13
G07090	C-TQTSY-1400	14	41	4	18	7	0.038	80	6, 10, 16
G05462	C-TQTSY-1600	16	41	4	18	7	0.038	80	8, 12, 18
G05461	C-TQTSY-1800	18	41	4	18	7	0.038	80	8, 14, 20
G05463	C-TQTSY-2000	20	41	4	18	7	0.038	80	8, 14, 22
G04220	C-TQTSY-2400	24	41	4	18	7	0.038	80	12, 18, 26
G04226	C-TQTSY-2800	28	41	4	18	7	0.038	80	14, 20, 30
G07242	C-TQTSY-3200	32	41	4	18	7	0.038	80	14, 22, 34






Fuhrman Pleural/Pneumopericardial Drainage Set – Seldinger, Pigtail

Used for evacuation of air from the pericardial sac or drainage of air or fluid from the pleural space.

Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- Sideports are positioned within the catheter pigtail and are designed to aid in drainage.
- Centimeter markings on the pigtail catheter guide insertion.
- The catheter material has radiopaque properties.
- Pigtail catheters have been shown to have a comparable efficacy to chest tubes in patients with pneumothorax.¹

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr
G55713	C-PPD-500-WCE-IMH*	5	15	6	18	5	0.028	40	4.5
G55714	C-PPD-600-WCE-IMH*	6	15	6	18	5	0.028	40	5.5
G55723	C-PPD-850-IMH	8.5	15	6	18	7	0.028	40	8
G55715	C-PPD-850-WCE-IMH*	8.5	15	6	18	5	0.028	40	8
G55716	C-PPD-1020-WCE-IMH*	10.2	15	6	18	5	0.028	40	9.5
G55727	C-PPD-1200-CHILDRENS-082694-IMH	12	15	6	18	7	0.028	40	11



^{1.} Kulvatunyou N, Vijayasekaran A, Hansen A, et al. Two-year experience of using pigtail catheters to treat traumatic pneumothorax: a changing trend. *J Trauma*. 2011;71(5):1104-1107.



Fuhrman Pleural/Pneumopericardial Drainage Tray – Seldinger, Pigtail

Used for evacuation of air from the pericardial sac or drainage of air or fluid from the pleural space.


Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- Sideports are positioned within the catheter pigtail and are designed to aid in drainage.
- Centimeter markings on the pigtail catheter guide insertion.
- The catheter material has radiopaque properties.
- Pigtail catheters have been shown to have a comparable efficacy to chest tubes in patients with pneumothorax.¹

1. Kulvatunyou N, Vijayasekaran A, Hansen A, et al. Two-year experience of using pigtail catheters to treat traumatic pneumothorax: a changing trend. *J Trauma*. 2011;71(5):1104-1107.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr
G55722	C-PPDY-700-022394-IMH	7	15	6	18	7	0.028	40	6.5
G55724	C-PPDY-850-IMH	8.5	15	6	18	7	0.028	40	8
G55726	C-PPDY-1020-OH-011795-IMH	10.2	25	6	18	7	0.028	60	8, 10

 $Some\ products\ or\ part\ numbers\ may\ not\ be\ available\ in\ all\ markets.\ Contact\ your\ local\ Cook\ Medical\ sales\ representative\ or\ Customer\ Support\ \&\ Distribution\ for\ details.$



Additional components

- 5 mL of lidocaine
- 22 and 25 gage Monoject needles
- 3 and 6 mL Monoject syringes
- #11 thumb scalpel
- 30 x 30 inch CSR wrap

- 24 x 36 inch fenestrated drape
- ChloraPrep
- 4 x 4 inch gauze sponges (pack of 10)
- Braided black silk suture
- Catheter retention disk
- Ty-Rap cable tie
- Mosquito Serrated Hemostat



Lock Pericardiocentesis Catheter Set and Tray – Pigtail Catheter

Used to remove fluid from the pericardial sac.


Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The 40 cm long catheter is designed to cater to larger patient anatomies.
- The radiopaque catheter material enhances x-ray visualization.

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Dilator Fr
Set									
G02433	C-PCS-830-LOCK	8.3	40	6	18	7, 12	0.038	70	8
Tray									
G04507	C-PCSY-830-LOCK	8.3	40	6	18	7, 12	0.038	70	8

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Support & Distribution for details.

Set components





Peritoneal Lavage Set and Tray

Used to infuse or withdraw rinsing solutions in the peritoneal cavity.

Features and benefits


- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The catheter is equipped with 90 sideports.
- The radiopaque catheter material enhances x-ray visualization.


Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm	Wire Guide Diameter inch	Wire Guide Length cm
Set								
G05267	C-PLS-901	9	20	90	18	7	0.038	40
Tray								
G05268	C-PLSY-901	9	20	90	18	7	0.038	40

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Support & Distribution for details.

Set components







PTFE-COATED FIXED CORE WIRE GUIDE

Additional tray components

- 5 mL of lidocaine
- Braided black silk suture
- 22 and 25 gage Monoject needles
- 3 and 6 mL Monoject syringes
- #11 thumb scapel
- 4 x 4 inch gauze sponges (pack of 10)
- Needle holder cup
- ChloraPrep
- 30 x 30 inch CSR wrap
- 18 x 24 inch fenestrated drape



Simple Pneumothorax Aspiration Accessory Set

Used in conjunction with a pneumothorax aspiration catheter if incomplete expansion of the lung persists after mechanical aspiration of air.


Features and benefits

- The set comes with the Cook Chest Drain Valve, which is designed to help remove air from the pleural cavity.
- The Cook Chest Drain Valve can be used for the transportation of patients who need to have pleural drainage catheters in place for long durations, or for outpatient procedures.

Order Number	Reference Part Number	Connecting Tube Connecting Tube Length Fr cm		
Set				
G05297	C-CASP-A-FORD	14 30		
Cook Chest Drain Valve Only				
G36370	C-CDV-1	-		

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Support & Distribution for details.

VINYL CONNECTING TUBE





Plastic Multipurpose Tubing Adapter

Used for connection to other medical devices such as catheters.

Features and benefits

• The adapter is a male Luer Lock to a universal taper.

Order Number	Reference Part Number	Compatible Instrumentation Fr	
G05302	C-PMLLA-VTA-L	20-40	

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Support & Distribution for details.



Vinyl Connecting Tube

Used for the connection of an external drainage catheter to a drainage pouch.

Order Number	Reference Part Number	Connecting Tube Fr	Connecting Tube Length cm	
Stopcock				
G02898	C-CTU14.0-30-ST	14	30	
G02327	CTU14.0-40-ST	14	40	
Luer Lock				
G02791	CTU14.0-30	14	30	



Thal-Quick Chest Tube Adapter

Used to provide separate sampling or infusion port within Thal-Quick Chest Tubes.

Features and benefits

• This accessory is compatible with the Thal-Quick Chest Tube Set and Tray (see pages 10 and 11).

Order Number	Reference Part Number			
G07467	C-TQF-28.0			

Customer Service

EU Website: cookmedical.eu EDI: cookmedical.eu/edi

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com Belgium: +32 27001702, be.orders@cookmedical.com Denmark: +45 38487607, da.orders@cookmedical.com Finland: +358 972519996, fi.orders@cookmedical.com France: +33 171230269, fr.orders@cookmedical.com Germany: +49 6950072804, de.orders@cookmedical.com Hungary: +36 17779199, hu.orders@cookmedical.com

Iceland: +354 8007615, IS.orders@cookmedical.com Ireland: +353 61239252, ie.orders@cookmedical.com Italy: +39 0269682853, it.orders@cookmedical.com Netherlands: +31 202013367, nl.orders@cookmedical.com Poland: +48 223060159, pl.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customer support @cook medical.com

Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com





CC, URO-8.5X11