

BriteGrip™

Lightweight Disposable Laryngoscope Handle



The Performance & Power of a Reusable Laryngoscope Handle without the Weight, Contamination Hazard or Hassle


Weights Only 3oz

The lightweight Hartwell Medical, single-use disposable laryngoscope handle provides the performance and power of a reusable. This innovative handle offers numerous advantages to substantially reduce the potential for cross-contamination known to exist with reusable devices. Plus, you'll save valuable time and hassle never needing to sterilize and inspect your handles again.

Made from high-quality materials and a sturdy construction insures reliability in the most adverse conditions. An ideal companion product to BRITVIEW® and GRANDVIEW® laryngoscope blades.

Features and Benefits

- **Incredibly Lightweight**
- **Batteries Included, Ready To Go!**
- **Single-Use Design & Disposable**
- **Reduces Potential Cross-Contamination**
- **Sturdy Construction**
- **Powerful & Reliable and Recyclable**
- **Eliminates Equipment Sterilization**
- **Compatible With Standard Blades That Conform To ASTM (F965) And ISO (7376-1)**



www.HartwellMedical.com



HARTWELL MEDICAL

6354 Corte del Abeto, Suite F
Carlsbad, CA 92011
(800) 633-5900
(760) 438-5500
Fax: (760) 438-2783