

UCapIt™ CAP 3 and CAP 5 Controlled Access Pharmaceutical Dispensers

06.23



The CAP 3 and CAP 5 Controlled Access Pharmaceutical Dispensers are customizable inventory management tools that help address major supply room challenges.

CAP dispensers release items to EMS personnel based on customizable ID protocols. They can also prompt users for station or dispatch specific data. All supply withdrawals are cataloged, time stamped and regulated by tech clearance/training level.

EMS management can monitor and restock inventory at any time. CAP machines aggregate dispense history, create pick lists for restocking and generate individual call reports based on products restocked. Managers can choose to receive automatic email or text alerts regarding low stock, temperature variance, usage spikes, expiration dates or power loss.



FEATURES

- Provides customizable dispensing and reporting options for individual agency needs
- Accepts a variety of ID protocols including individual pin numbers, proximity cards, magnetic stripe cards, bar codes, biometrics or a combination of these methods
- Monitors and tracks access to controlled substances, narcotics and pharmaceuticals
- Helps ensure that items with expiration dates are used in a timely manner

AVAILABLE ADD-ONS

- Refrigeration
- Lock down security package
- Extension machine (no controller)

UCapIt™ CAP 3 and CAP 5 Controlled Access Pharmaceutical Dispensers

06.23



PRODUCT SPECIFICATIONS

Included features:

- Vibrant 7" Color Touchscreen
- iVend® Delivery Sensor System
- Parrot Door (enables re-use of protective bags)
- LED lighting
- ADA Compliant

Optional Features:

- Proximity badge reader
- Biometric reader*

*Requires a registration station



CAP 3 SPECIFICATIONS

- 23 standard selections
- Holds 384 standard items
- Electrical requirements: 115 VAC/60Hz 1.2 AMPS
- Weight: 445 lbs
- 72" x 29" x 34"



CAP 5 SPECIFICATIONS

- 40 standard selections
- Holds 630 standard items
- Electrical requirements: 115 VAC/60Hz 1.2 AMPS
- Weight: 628 lbs
- 72" x 41" x 34"

ITEM #	DESCRIPTION
4800-001H	CAP 3 Standard
4310-75714	Biometric Reader
4310-75716	Biometric Registration Station
4800-31184	Prox Card Reader
4800-43110	iQ Technology 12 Month Subscription

ITEM #	DESCRIPTION
4800-0001H	CAP 5 Standard
4800-0002S	CAP 5 Extension Machine
4800-0-03-R	CAP 5 Refrigerated Machine - Bottom Half
4310-75714	Biometric Reader
4310-75716	Biometric Registration Station
4800-31184	Prox Card Reader
4800-43110	iQ Technology 12 Month Subscription

UCapIt™ CAP Adjustable Locker

06.23



The UCapIt™ CAP Adjustable Locker is an inventory management tool that helps address major supply room challenges. It can be used alone or connected to the CAP 5 or CAP 3 Controlled Access Pharmaceutical Dispenser.

The locker is ideal for dispensing large, bulky items or items that need to be checked out and checked back in by EMS personnel. Items can be checked out based on customizable ID protocols. All supply withdrawals are catalogued, time stamped and regulated by tech clearance/training levels. Lockers are able to be quickly and easily configured on-site to accommodate various product mixes with up to 20 compartments.

SPECIFICATIONS

- Up to 20 Compartments
- All compartments are 10 ¹⁵/₁₆" W x 22 ¹⁵/₁₆" D; Height is variable
- Electrical req: 115 VAC/60Hz, 1.2 AMPS
- 73" H x 30" W x 26" D
- 327 lbs*

*Weight varies with configuration



ITEM #	DESCRIPTION
4800-01H	Stand-Alone Adjustable Locker w/Controller, BTM Black
4800-02S	Satellite Adjustable Locker, BTM Black