

WHY CHOOSE THE V+PRO VENTILATOR FOR EMS?

3-in-1 Device

Includes Invasive, Non-Invasive, and integrate HFNC (ranging from 4-60 L/min), right out of the box.

9 Modes of Ventilation

- | | |
|--|--|
| 1. Spontaneous | |
| 2. Bi-Level | 7. Vol. Targeted-PS |
| <i>with functionality similar to S/T, Timed, and BiPAP ventilation</i> | <i>with functionality similar to AVAPS® and iVAPS®</i> |
| 3. Assist/Control-Pressure | 8. Vol. Targeted-PC |
| 4. Assist/Control-Volume | <i>with functionality similar to PRVC</i> |
| 5. SIMV-Pressure <i>including CPAP</i> | 9. Vol. Targeted-SIMV |
| 6. SIMV-Volume | <i>with functionality similar to SIMV+PRVC</i> |

Robust Leak Compensation

The Leak+ feature allows VOCSN to deliver High Flow ventilation as well as compensate for a leak up to 175 L/min.

Hospital to Home

Approved for patient transport both inside and outside the hospital.

Wide Range for Patient Care

Suitable for pediatric through adult patient population weighing at least 5 kg.

Long Battery Life

Contains 2 removable 4-hour batteries and 1 internal 1-hour battery for a total battery life of 9 hours.

High Standard of Care

Meets rigorous ISO 80601-2-12 international standard for critical-care ventilators.

DISPOSABLE PARTS

Adult Passive Vent Circuit
PRT-00946-000

Adult Active Vent Circuit
PRT-00950-001

Ped. Active Vent Circuit
PRT-01011-000

Ped. Passive Vent Circuit
PRT-01015-000

External Bacterial Filter (10-Pack)
PRT-00850-000

