



Curaplex® Safety Lancets are a top choice for capillary blood collection fingerstick users: ease of use, comfort, safety and reliability!

- Easy to use
- Minimizes pain
- Safe, single-use design
- Proven reliability
- 21, 25, 28 Gauge needles
- Eco-friendly design
- 2.2 mm depth penetration

Ordering Information

2764-70721 - 21 Gauge, Green, 100/BX, 48BX/CS

2764-70725 - 25 Gauge, Blue, 100/BX, 48BX/CS

2764-70728 - 28 Gauge, Orange, 100/BX, 48BX/CS

For more information, please contact your Account Manager or call Customer Service at 800.533.0523



Bound Tree

800.533.0523 | www.boundtree.com

Curaplex® Safety Lancets



CROSS REFERENCE CHART

Brand (Vendor)	Item #	Gauge	Depth (mm)	Curaplex Item #	Gauge	Depth (mm)	Flow	Color
Medi Purpose (Surgilance)	E6255	21G	2.8	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Medi Purpose (Surgilance)	E6254	21G	2.2	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Medi Purpose (Surgilance)	782-SLN200BX	21G	1.8	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Emergency Medical Products Inc.	2764-12107	21G	1.8	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Medline Industries, Inc.	2764-26021	21G	2.2	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Owen Mumford Inc.	675-AT0702	21G	2.4	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Owen Mumford Inc.	675-AT0712	21G	3.0	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Owen Mumford Inc.	2764-10141	21G	2.0	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Owen Mumford Inc.	2764-10121	21G	2.0	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Med Plus Services Usa	2764-10080	21G	2.4	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Htl Strefa, Inc.	970121	21G	1.8	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Dynarex Corporation	7112	21G	2.2	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Medipurpose	SLN300	21G	2.8	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Owen Mumford, Inc.	AT1012	21G	2.0	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Medipurpose	SLN240	22G	2.2	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Owen Mumford, Inc.	AT1002	23G	1.8	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Medipurpose	SLN200	23G	1.8	2764-70721	21G	2.2	Medium - 20-40 µL	Green
Htl Strefa, Inc.	970125	25G	1.4	2764-70725	25G	2.2	Medium - 20-40 µL	Blue
Dynarex Corporation	7114	26G	1.8	2764-70725	25G	2.2	Medium - 20-40 µL	Blue
Dynarex Corporation	1880-71426	26G	1.8	2764-70725	25G	2.2	Medium - 20-40 µL	Blue
Owen Mumford Inc.	675-AT1042	28G	1.8	2764-70728	28G	2.2	Low - 10-20 µL	Orange
Owen Mumford Inc.	55845AT0714	28G	1.8	2764-70728	28G	2.2	Low - 10-20 µL	Orange
Owen Mumford Inc.	2764-10441	28G	1.8	2764-70728	28G	2.2	Low - 10-20 µL	Orange
Medipurpose	SLN100	26G	1.0	2764-70728	28G	2.2	Low - 10-20 µL	Orange
Dynarex Corporation	7115	28G	1.8	2764-70728	28G	2.2	Low - 10-20 µL	Orange
Owen Mumford, Inc.	AT1042	28G	1.8	2764-70728	28G	2.2	Low - 10-20 µL	Orange
Cardinal Health Aurora	3007788	28G	-	2764-70728	28G	2.2	Low - 10-20 µL	Orange
Htl Strefa, Inc.	970128	28G	1.6	2764-70728	28G	2.2	Low - 10-20 µL	Orange
Owen Mumford, Inc.	AT1022	30G	1.5	2764-70728	28G	2.2	Low - 10-20 µL	Orange

For more information, please contact your Account Manager or call Customer Service at 800.533.0523



Bound Tree

800.533.0523 | www.boundtree.com