

# Curaplex® Emesis Bags and Dispenser



**our brand.  
your advantage.**



**Cost-effective options for medical providers in a variety of roles. The convenient, mountable dispenser provides a secure holding for the emesis bags.**

## Emesis Bag features:

- Offered in a transparent HDPE (High Density) material or LDPE material, which is thicker, softer, and more flexible
- Sickness bags are designed to reduce contamination and smell
- Bag has a patient convenient easy lock style designed for additional securement
- Rim made of durable plastic
- Tear resistant bag
- Graduated markings to measure oz. and cc/ml
- Not made with natural rubber latex
- Single use

## Emesis Bag Dispenser features:

- No assembly required for the emesis bag dispenser
- Easy to mount to wall with included hardware
- Made of durable polyethylene and rubber
- Top-lid covers emesis bags, aiding in infection control
- Securely holds up to 25 folded emesis bags

ITEM #	DESCRIPTION	QTY
<b>1071-11106</b>	Dispenser, Holds 25	1/ea, 6/cs
<b>1071-10204</b>	1000cc, Low Density	25/pk, 6 pk/cs
<b>1071-10202</b>	1000cc, High Density	25/pk, 6 pk/cs