

Curaplex® IV Administration & Extension Sets



Curaplex needle-free administration and extension sets come in a variety of configurations and lengths to meet your needs for patient procedures.

FEATURES:

- Individually packaged
- Sterile
- Not made with natural rubber latex
- DEHP free
- Non-pyrogenic



Item # 1714-31081

Priming Volume: 2.0mL

50 ea/cs

Curaplex® 8in Extension Set, Removable Sure-Lok, NeedleFree Connect, Pinch Clamp, Rotating Male LL




Item # 1712-10830

Priming Volume: 16mL

50 ea/cs

Curaplex® I.V. Admin Set, 10 Drop, 83in, Pre-Pierced Y-Site, Sure-Lok Needle-Free Y-Site, Rotating Male Luer Lock



Curaplex® IV Administration & Extension Sets




Item # 1712-10835

Priming Volume: 18mL

50 ea/cs

Curaplex® I.V. Admin Set, 10 Drop, 89in, Roll Clamp, PP Y-Site, Sure-Lok Needle-free Y-Site, Attached Extension Set



Item # 1712-15830

Priming Volume: 18mL

50 ea/cs

Curaplex® I.V. Admin Set, 15 Drop, 89in, Roll Clamp, PP Y-site, Sure-Lok Needle-Free, Attached Extension Set




Item # 1712-60830

Priming Volume: 15mL

50 ea/cs

Curaplex® I.V. Admin Set, 60 Drop, 83in, PP Y-Site, Sure-Lok Needle-free Y-Site, Rotat Male Luer Lock



Item # 1812-02001

Priming Volume: 0.04mL

100 ea/cs

Curaplex® Needle-Free Male Adapter Plug

