

# The Non-Contact Infrared Thermometer



## IR 300

This is an FDA 510k cleared device

**Mode Indicator**

**Easy To Understand Digital Screen**

The screen background will illuminate to indicate no fever, slight fever, or high fever.

<p>95.9 °F - 99.2 °F No Fever</p>	<p>99.3 °F - 100.5 °F Slight Fever</p>	<p>100.6 °F - 109.2 °F High Fever</p>
---------------------------------------	--	---



**Set Button**  
Choose between audible or silent mode.

**Memo Button**  
View up to 50 previous readings.

**Mode Button**  
Two temperature modes are available: Body and Surface.

**On/Scan Button**  
Easy one handed operation.

## SPECIFICATIONS:

- **Measurement Time:** ≤ 2 seconds
- **Measuring Distance:** 0.4" - 2.0"
- **Body Temperature Range:** 95.0 °F - 107.6 °F
- **Ambient Temperature:** 59 °F - 95 °F
- **Accuracy:** ± 0.4 °F
- **Automatic Power-Off:** ≤ 18 seconds
- **Batteries:** 2 x AAA
- **Product Longevity:** 13 years
- **Product Weight:** 3 ounces (without batteries)
- **Product Size:** L 5.43" x W 3.74" x H 1.47"

ITEM #	DESCRIPTION	UOM
2731-31002	Non-Contact Infrared Thermometer, IR 300	1/EA, 40 EA/CS