

TRUST WHAT'S IN YOUR HANDS

Nonin CO-Pilot™

Wireless Handheld Multi-Parameter System (H500)



Pulse Oximetry



Pulse Rate



Carboxyhemoglobin



Methemoglobin



Cerebral or Tissue
Oximetry

**Emergencies can happen anywhere.
And unreliable measurements leave everyone exposed.**

You need to be able to depend on the tools you use so you can focus on the lives that depend on you. Introducing the Nonin CO-Pilot™ Multi-Parameter Wireless Handheld System to reliably measure five parameters including SpO2, COHb, and rSO2 in even the most challenging conditions. **So you can trust what's in your hands.**






A First-of-its-Kind Wireless, Multi-Parameter Monitoring System

Designed For Demanding EMS Environments

Engineered to withstand challenging conditions faced in the field, the Nonin CO-Pilot™ System is a flexible, five-parameter system capable of measuring SpO₂, Pulse Rate, COHb, MetHb and rSO₂ all on one wireless handheld device.

Responsive Readings

-  Accepts Nonin SpO₂, rSO₂ and multi-parameter sensors
-  Indicated for SpO₂ accuracy from 80% to 100% with CO up to 15%
-  COHb accuracy with patients in normoxia and hypoxia

Reliability in the Field

-  Lightweight, handheld design is easy to carry and simple to use
-  Nonin multi-parameter and rSO₂ sensors deliver dependable performance
-  Rugged design meets strict IEC 60601-1-12 EMS standards

Features Designed to Give Accuracy and Reliability

Proven durable sensors can monitor multiple parameters at once

Pleth waveform displays simultaneously with the parameters measured using the Multi-Parameter Sensor (8330AA)

Wireless display and signal processor are pre-paired; no additional steps needed

System can also monitor rSO₂ during resuscitation or ECMO

Full battery charge provides approximately 10 hours of runtime



Nonin Medical, Inc.

13700 1st Avenue North
Plymouth, MN • 55441-5443 • U.S.A.

Tel: +1.763.553.9968 1.800.356.8874
Fax: +1.763.577.5521
Email: info@nonin.com

Nonin Medical B.V.

Doctor Paul Janssenweg 150
5026 RH Tilburg • Netherlands

Tel: +31 (0)13-79 99 040
Fax: +31 (0)13-79 99 042
Email: infointl@nonin.com



nonin.com

CE 0123



Rx Only