


CuraSafe™ Comfort BC Push Button Safety IV Catheter/Blood Control

VenaPUR® (polyurethane) Material designed to remain kink-free

product specifications

- Incorporated blood control
- VenaPUR® catheter material
- Push button style needle retraction
- Ground stainless-steel needle
- Straight hub
- Disposable
- Non-pyrogenic
- Not made with natural rubber latex
- Not made with DEHP
- ISO 10555-5 compliant color-coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant



ordering information

Color	Item #	Gauge and Length	Catheter ID (mm)	Catheter OD (mm)	Gravity Flow Rate (mL/min)	Qty
Orange	1613-86414	14G x 1.75"	1.6mm	2.1mm	297mL/min	50/bx, 200/case
Grey	1613-86412	16G x 1.16"	1.38mm	1.8mm	220mL/min	50/bx, 200/case
Green	1613-64162	18G x 1.88"	0.95mm	1.3mm	94mL/min	50/bx, 200/case
Green	1613-86418	18G x 1.16"	0.95mm	1.3mm	94mL/min	50/bx, 200/case
Pink	1613-64202	20G x 1.88"	0.78mm	1.1mm	61mL/min	50/bx, 200/case
Pink	1613-86420	20G x 1.16"	0.78mm	1.1mm	61mL/min	50/bx, 200/case
Pink	1613-64203	20G x 1"	0.78mm	1.1mm	61mL/min	50/bx, 200/case
Blue	1613-86422	22G x 1"	0.62mm	0.9mm	36mL/min	50/bx, 200/case
Yellow	1613-86424	24G x .75"	0.49mm	0.7mm	20mL/min	50/bx, 200/case

