


Endotracheal Tubes without Stylette



These endotracheal tubes are constructed of kink-resistant medical-grade PVC and feature a radio-opaque black line. Available with or without a cuff. The pilot balloon on the cuffed versions allows for a tactile sense of cuff volume. Latex free and packaged sterile.

Ordering Information

Uncuffed

2113-20320	2.0 mm
2113-20325	2.5 mm
2113-20330	3.0 mm
2113-20335	3.5 mm
2113-20340	4.0 mm
2113-20345	4.5 mm
2113-20350	5.0 mm
2113-20355	5.5 mm
2113-20360	6.0 mm
2113-20365	6.5 mm

Cuffed

2113-20230	3.0 mm	2113-20270	7.0 mm
2113-20235	3.5 mm	2113-20275	7.5 mm
2113-20240	4.0 mm	2113-20280	8.0 mm
2113-20245	4.5 mm	2113-20285	8.5 mm
2113-20250	5.0 mm	2113-20290	9.0 mm
2113-20255	5.5 mm	2113-20295	9.5 mm
2113-20260	6.0 mm	2113-20200	10.0 mm
2113-20265	6.5 mm		

Product Specifications

- Standard 15 mm proximal connection with Murphy tip
- Sterile packaging for single-patient use
- Cuffs sized proportionally to the tube size
- Latex free

For more information, please contact your Account Manager or call Customer Service at 800.533.0523



Bound Tree

800.533.0523 | www.boundtree.com