


# Endotracheal Tubes with Stylette



This set includes a quality endotracheal tube with an inserted stylette, saving time and storage space. The stylette has a low-friction coating to facilitate smooth withdrawal. It is designed with a soft distal tip to reduce the risk of trauma. The tubes are available with or without a cuff. The pilot balloon allows for tactile sense of cuff volume.

## Product Specifications

- Standard 15 mm proximal connection with Murphy tip
- Sterile packaging for single-patient use
- Cuffs sized proportionally to the tube size
- Latex free

## Ordering Information

Uncuffed		Cuffed			
2113-10325	2.5 mm	2113-10225	2.5 mm	2113-10265	6.5 mm
2113-10330	3.0 mm	2113-10230	3.0 mm	2113-10270	7.0 mm
2113-10335	3.5 mm	2113-10235	3.5 mm	2113-10275	7.5 mm
2113-10340	4.0 mm	2113-10240	4.0 mm	2113-10280	8.0 mm
2113-10345	4.5 mm	2113-10245	4.5 mm	2113-10285	8.5 mm
2113-10350	5.0 mm	2113-10250	5.0 mm	2113-10290	9.0 mm
		2113-10255	5.5 mm	2113-10200	10.0 mm
		2113-10260	6.0 mm		

For more information, please contact your Account Manager or call Customer Service at 800.533.0523



**Bound Tree**

800.533.0523 | [www.boundtree.com](http://www.boundtree.com)