



DICK MEDICAL SUPPLY
1946 CAMARO AVENUE
COLUMBUS, OHIO 43207
614-444-2300 | WWW.DICKMEDICAL.COM

Technical Data Sheet

SureShape™ Moldable Splint 4.25"x9"



Item #64009

Contact Information:

Distributor: Dick Medical Supply, LLC.
Address: 1946 Camaro Ave
Columbus, Ohio 43207 USA
Phone: (614) 444-2300
Katrina Copeland
E-Mail: kcopeland@dms4ems.com
Revision: 02152022



DICK MEDICAL SUPPLY
1946 CAMARO AVENUE
COLUMBUS, OHIO 43207
614-444-2300 | WWW.DICKMEDICAL.COM

1. Product Overview

The SureShape™ moldable splint is a highly versatile, lightweight medical device used for immobilizing bone and soft tissue injuries. Often referred to as a "pocket cast," it is staple in emergency medicine, wilderness survival, and first aid kits due to its unique balance of flexibility and strength.

- **Construction:** A thin core of soft aluminum alloy sandwiched between two layers of closed-cell foam.
- **Radiolucent:** It does not need to be removed for X-rays.
- **Waterproof:** The foam exterior is non-absorbent and easy to clean.
- **Compact:** It can be rolled or folded for storage in small trauma kits.

2. Specifications

Product Name:	SureShape™ Moldable Splint, 4.25"x9", Each
Color:	Orange
Size:	Universal
Certification:	CE, FDA
Reusable:	Yes
Country of Origin:	USA
Application:	EMS
Item#:	64009
Pack Size:	EACH

Disclaimer: Dick Medical Supply, LLC believes that the information in this Technical Data Sheet is accurate and represents the best and most current information available to us. Dick Medical Supply makes no representations or warranties either express or implied, regarding the suitability of the material for any purpose or the accuracy of the information contained within this document. Accordingly, Dick Medical Supply will not be responsible for damages resulting from use of or reliance upon this information.