

■ **SIMBODIES® offer the most lifelike medical manikins available for medical simulation.** While SIMBODIES manikins provide an immersive simulation experience, they may have limitations based on the specific technical medical intervention being practiced.

Intended Medical Procedures	Adult				Pediatric			Infant		Lifeless
	Pro	EMS-Man	EMS-TT	Standard	Advanced (10 y/o)	Standard (10 y/o)	Advanced (5 y/o)	Advanced	Pre-term (28 wk)	Burnt & Blasted
Catastrophic / Life-Threatening Hemorrhage										
Tourniquet application	✓	✓		✓	✓	✓	✓			
Wound packing	✓	✓		✓*	✓*					
Junctional bleeding	✓			✓*	✓*					
Airway Management										
Surgical airway (cricothyroidotomy)	✓		✓							
Airway examination	✓	✓	✓		✓		✓	✓	✓	
Head-tilt / chin-lift / jaw-thrust	✓	✓	✓		✓		✓	✓	✓	
Airway suction	✓	✓	✓		✓		✓	✓	✓	
Postural / nasopharyngeal / oropharyngeal airway management	✓**	✓**	✓**		✓**		✓**		✓**	
Supraglottic airway device insertion	✓	✓	✓		✓		✓	✓		
Laryngoscopy / tracheal intubation	✓	✓	✓		✓		✓	✓	✓	
Maxillo-facial packing	✓	✓	✓		✓		✓	✓		
C-Spine Control										
Spinal collar application	✓	✓	✓	✓	✓	✓	✓	✓		
Manual inline stabilization application	✓	✓	✓	✓	✓	✓	✓	✓		
Breathing Assessment										
Artificial ventilation – BVM	✓	✓	✓		✓		✓	✓		
Artificial ventilation – intubated	✓	✓	✓		✓		✓	✓		
Chest seal application	✓	✓	✓	✓	✓	✓	✓	✓		
Finger thoracostomy site	✓		✓							
Needle thoracentesis	✓	✓	✓		✓		✓	✓		
Circulation Assessment										
Dressing application (various)	✓	✓		✓	✓	✓	✓	✓		
Pelvic binder application	✓	✓		✓	✓	✓	✓	✓		
Traction device application – femoral	✓	✓		✓	✓	✓	✓	✓		
Intraosseous access – tibia	✓	✓			✓		✓	✓		
Intraosseous access – humeral head	✓	✓	✓		✓		✓			
Chest compressions – manual	✓	✓	✓	✓	✓	✓	✓	✓	✓	
Chest compressions– automated	✓**	✓**	✓**		✓**		✓**	✓**		
AED application / training defibrillator pads	✓	✓	✓	✓	✓	✓	✓	✓		✓

* Bespoke option

** For most current medical devices