



eitan
medical
sapphire™

Sapphire™ FastTest™ PM

Your solution for independent
onsite infusion pump
maintenance



Want to streamline and automate your facility's maintenance process for infusion pumps and avoid relying on outside service centers?

FasTest PM is the solution

FasTest PM is an on-site maintenance servicing tool for the Sapphire Infusion Pump Platform, which enables you to:



Quickly and easily conduct full maintenance verification tests, including occlusion detection, air detection, and pump accuracy



Reduce maintenance costs and equipment downtime, which puts you in control



Receive an immediate certification report helping maintain compliance





FasTest PM gives you independence and control of the maintenance process

Standardize pump testing with a dedicated software tool. FasTest PM offers you the option to do annual PM tests between treating patients, keeping more pumps available, and verified for patient use.

FasTest PM makes pump maintenance easier than ever

Utilize this easy and intuitive solution, which provides step-by-step, guided navigation.

FasTest PM brings efficiency home

Save costs, time, and improve efficiency by moving your annual PM tests in-house. No more chasing down pumps. Within minutes, you can have your manufacturer certification, ensuring your infusion devices are verified.

FasTest PM helps you streamline your maintenance processes

Deploy an efficient maintenance system that is fast, cost-effective and sustainable. FasTest is reusable and offers multiple applications, giving your team full independence to conduct maintenance on your Sapphire fleet.



Technical Features

No. of tests per kit

10

Sequence of tests

Visual inspection, Safety features testing, Calibration verification

PM run time

≈ 12 min

Maximum simultaneous PM processes

8

Calibration verification process

- Air Detection
- Flow Accuracy
- Occlusion Sensing

Kit content

- Administration set - 1800 mm tubing, with: spike, syringe, back check valve, Y-connector, T-connector, administration cassette, slide clamp, sleeve, and fluid bag cap
- 500 mL fluid bag

FasTest PM certificate

Automatically generated certificate

- Device SW version
- Certification kit S/N
- Test date
- Technician name
- Visual inspection & safety features test results
- Air Detection, accuracy & occlusion test results

16024-003-0007 Rev 07 | 09.2025

Contact Information:

service@eitanmedical.com

877.541.9944



www.eitanmedical.com