

SAPPHIRE
INFUSION PLATFORM
US CATALOG 2020



SAPPHIRE
by Eitan Group

CONTENTS

●	<u>INFUSION KITS</u>	4
●	<u>ACCESSORIES</u>	8
●	<u>ADMINISTRATION SETS</u>	12
●	<u>PC APPLICATIONS</u>	15

Sapphire™ Multi-therapy infusion pump kit | SW Rev 15.00.0



Includes:

- Sapphire™ Multi-therapy pump
- Power supply
- Mini cradle
- Mini cradle knob key
- User manual

Power supply type:

- Fixed head power supply
Kit P/N 17000-028-0072
- Integrated power supply & AC Power Cord - Type B
Kit P/N 17000-028-0077
- Interchangeable plug power supply
Kit P/N 17000-028-0078

Sapphire™ Epidural infusion pump kit | SW Rev 15.00.0



Includes:

- Sapphire™ Epidural pump
- Power supply
- Bolus handle
- Mini cradle with splitter
- Mini cradle knob key
- User manual

Kit P/N 17000-031-0035

Sapphire™ Multi-therapy infusion pump kit | SW Rev 13.23.2

Sapphire™ Epidural infusion pump kit | SW Rev 13.23.2



Includes:

- Sapphire™ Multi-therapy pump
- Power supply
- Mini cradle
- Mini cradle knob key
- User manual

Kit P/N 17000-028-0061

Includes:

- Sapphire™ Epidural pump
- Power supply
- Bolus handle
- Mini cradle with splitter
- Mini cradle knob key
- User manual

Kit P/N 17000-028-0027

ACCESSORIES

Sapphire™ Mini Cradle

Sapphire™ Mini Cradle

With Splitter

P/N 17000-026-0007

Without Splitter

P/N 15042-000-0010



Sapphire™ Multi-Pump Mounting System



P/N 17000-020-0008

Sapphire™ Power Supply



P/N 15072-000-0007

Sapphire™ External Battery Pack



P/N 15079-111-0001

Sapphire™ Splitter



P/N 15043-000-0005

Sapphire™ Integrated Power Supply



P/N 15072-000-0010

Sapphire™ AC Power Cord - Type B



P/N 15070-112-0005

Sapphire™ Power Cord - compatible with IPS



P/N 15171-112-0002

ACCESSORIES

Sapphire™
Infusion Pouch 500 mL



P/N 05080-605-0001

Sapphire™
Travel Case



P/N 15020-010-0002

Sapphire™
Homecare Large Backpack



P/N 15088-000-0001

Sapphire™
PCA lockbox 250mL



P/N 15063-000-0001

Sapphire™
Lockbox 250mL strap and pouch



P/N 15063-000-0004

Sapphire™
PCA lockbox 100mL







P/N 15061-131-0004

Sapphire™
Bolus handle







P/N 15071-000-0005






ADMINISTRATION SETS | STANDARD SETS

NAME	P/N	DESCRIPTION	LENGTH	PRIME VOLUME	UNITS/BOX
AP 403-01	12003-000-0013	Microbore infusion set with non-vented spike	290 cm 114 in	~6 mL	30
					
AP 409-01	12005-000-0003	Macrobore primary infusion set with vented/non-vented drip chamber, back check valve, roller clamp and 2 needleless Y-site	277 cm 109 in	~23 mL	30
					
AP 423-01	12003-000-0023	Microbore infusion set with non-vented spike and back check valve	292 cm 115 in	~7 mL	30
					
AP 433-01	12000-000-0044	Macrobore infusion set with vented/non-vented drip chamber, roller clamp and needleless Y-site	286 cm 113 in	~22 mL	20
					

ADMINISTRATION SETS | FILTER SETS

NAME	P/N	DESCRIPTION	LENGTH	PRIME VOLUME	UNITS/BOX
AP 204-01	12003-000-0070	Microbore infusion set with non-vented spike, 0.2µ filter and back check valve	295 cm 116 in	~8 mL	30
					
AP 206-01	12000-000-0022	Macrobore infusion set with non-vented spike, 1.2µ filter	295 cm 116 in	~23 mL	20
					
AP 214-01	12004-000-0002	Microbore yellow marked infusion set with vented/non-vented spike and 0.2µ filter	295 cm 116 in	~8 mL	30
					
AP 240-01	12003-000-0069	Microbore infusion set with non-vented spike, 1.2µ filter and back check valve	295 cm 116 in	~8 mL	20
					

ADMINISTRATION SETS | SPECIALTY SETS

NAME	P/N	DESCRIPTION	LENGTH	PRIME VOLUME	UNITS/BOX
AP 133-01	12006-000-0010	Macrobore blood infusion set with 2 non-vented spikes, 200µ filter, 3 roller clamps and needleless Y-site	292 cm 115 in	30 mL	20
					
AP 402-01	12004-000-0004	Microbore yellow marked infusion set with non-vented spike	290 cm 114 in	~6 mL	30
					
AP 410-01	12005-000-0007	Macrobore secondary infusion set with vented/non-vented drip chamber. Includes an infusion bag hanger	71 cm 28 in	~9 mL	30
					
AP 424-01	12003-000-0024	Microbore infusion set with non-vented spike, side arm with female luer lock and back check valve	294 cm 116 in	~7 mL	30
					
AP 504-01	12000-000-0045	Macrobore UV-Resistant infusion set with vented/non-vented drip chamber, roller clamp and needleless Y-site	287 cm 113 in	~26 mL	20
					

PC APPLICATIONS

The Sapphire modern, software based platform expands beyond the pump to the PC applications suite that complements the pumps;



PUMP LOADER

Enables upgrading or reloading software to the Sapphire pump via a wire connection.



DRUG LIBRARY

Enables to set configuration settings, drug specific limit and more. When complete, the library can be uploaded to multiple Sapphire pumps.



FASTEST PM

(SW + disposable kit [16024-000-0002])
 FasTest PM is Preventative Maintenance Servicing Tool for Sapphire Infusion Pump platform. FasTest PM allows end-users to quickly and economically conduct full maintenance tests, including occlusion detection, air detection and pump accuracy checks. FasTest PM provides end-users with an immediate, downloadable authentic manufacturer certificate enhancing patient safety and regulatory compliance. The tool is designed to significantly improve maintenance costs and runtime.



ANNUAL CERTIFICATION KIT

(SW + disposable kit [15075-020-0001])
 Enables to run performance test (air detection, pressure detection and accuracy). Results are compared to manufacturing database per pump.







EVENT VIEWER










Enables downloading and viewing the event log of the Sapphire infusion pump on a computer.



ACCESSORIES | PC APPLICATIONS

PC APPLICATION	ITEM	P/N	PICTURE
FasTest PM	FasTest PM kit	16024-000-0002	
Annual Certification Test	Annual certification kit	15075-020-0001	
All	Communication cable	05020-110-0213	
All	RS232 to USB 4 port adapter	15077-000-0001	

SPARE PARTS | ACCESSORIES

SPARE PART FOR	ITEM	P/N	PICTURE
Sapphire pump	Li-Ion battery	05020-160-0001	
Mini cradle	Connection knob flange Key	05040-380-0019	
Mini cradle	Connection knob key	05040-380-0018	
Mini cradle	MC- AVP Connection Bolt Knob Tip	05040-380-0027	
Mini cradle	MC- AVP Press Plane	05040-380-0026	
PCA lockbox 100	Key	05070-190-0015	
PCA lockbox 250	Key	10010-122-0046	
PCA lockbox 250	Release clip (eliminates lockbox locking onto cradle)	15063-121-0029	
PCA Lockbox 500	Key	15062-132-0001	

REF 10050-003-0009 | Rev.04 / 04.2020

Distributed by

SOVEREIGN MEDICAL, INC.

**620H Valley Forge Rd
Hillsborough, NC 27278
919-644-1113
www.sovmed.com**

U.S. & Canada | Contact information:

65 Enterprise 3rd Floor
Aliso Viejo, CA 92656
Tel: 877-541-9944 | 949.303.7845
Fax: 949-288-5482
Email: customerservice@eitangroup.com
www.eitangroup.com



SAPPHIRE
by Eitan Group