

Ambu® ResCue Mask



Key Benefits

- Transparent dome allows rescuer to visually check for patient's lip color and vomitus
- Oxygen inlet (optional) facilitates increased oxygen delivery to patient
- Made of durable plastics for easy cleaning and long product life
- Available in hard case, soft case or plastic bag
- DEHP-free material

Ambu® ResCue Mask

The Ambu ResCue Mask is designed to provide mouth-to-mask ventilation.

As a physical barrier between the rescuer and the patient Ambu ResCue Mask is used to give protected rescue breathing for any patient (any age group) and may also be used with a resuscitator.

The Ambu ResCue Mask helps protect rescuers in an emergency situation. It prevents direct contact with patient's mouth, nose and face and helps to overcome hesitation to start resuscitation.

Specifications

Product	
Size	
Dimensions	13.8 x 11 x 4.8 cm
Weight hardcase	112 g
Weight softcase	67 g
Materials	
Mask	Polyvinylchloride
Softcase	Nylon
Hardcase	Polypropylene



Ambu ResCue Mask hardcase and softcase

Ordering Specifications

Item no.	Description
000 252 102	Ambu ResCue Mask with O ₂ inlet, elastic strap, and one-way valve, in hard case
000 252 104	Ambu ResCue Mask with O ₂ inlet, elastic strap, and one-way valve in nylon case
000 252 105	Ambu ResCue Mask does not include O ₂ inlet or elastic strap in soft case
000 252 120	Ambu ResCue Mask with O ₂ inlet and one-way valve in plastic bag
000 252 123	Ambu ResCue Mask adult and infant soft without O ₂ inlet in soft case
000 252 130	Ambu ResCue Mask American Red Cross in soft case
000 252 253	Ambu ResCue Mask one-way valve with filter (spare part)



Ambu ResCue Mask with adult and infant mask in a softcase

*ResCue Mask with O₂ inlet is Rx only in the USA.

Ambu USA
 6230 Old Dobbin Lane
 Columbia, MD 21045
 Tel. 800 262 8462
 Fax 800 262 8673
ambuUSA.com

US: Rx only

