Wi-Fi Simultaneous Image on Secondary Display

› Aitraq Max View 10.4" Display

› Phone, Tablet and PC



Additional wire connection to Airtraq MaxView Display and to most HDMI monitors

Airtraq +Plus

One VL fits **All** blade types & sizes



Blades made from RECYCLED Polycarbonate Your ECO Choice

Exclusively Distributed in the USA by

Mercury Medical® Superior Solutions

(800) 835-6633

www.mercurymed.com

Airtraq+Plus

All In One Eco

Video laryngoscope



VERSATILE

AFFORDABLE

SUSTAINABLE

Designed By:
Prodol Meditec S.A.

Intubation from Any Position in Any Scenario



Swivel



Flip





Image Rotates 180°

For Routine & Difficult Airways



Perfectly balanced

Full hermetic blades with integrated handle



3.5" touch screen. User friendly navigation



Video Laryngoscope into Charging Dock

Specifications

Reference	LDA600
Description	Video Laryngoscope
Screen	3.5" IPS. Resolution 640 x 480
Touch Panel	Viewing angle = 85º
Video Sensor	CMOS. Resolution: 800 x 800
Battery Capacity	150 min. Rechargeable
Wi-Fi	Range: 5 m (approx.)
Internal Memory	16 Gb. 10 hours (approx.)
IP Protection	IP6x (IEC 60529)
Weight	195 grams
Auto-Recording	Yes. Selected by user
Playback	Yes. Full functionality
Hybrid Intubation	Yes. Manual control of VL light

Single Use ECO Blades

ISCC International Sustainability & Carbon Certification

REF.	DESCRIPTION	SIZE
LDA6M4	MacIntosh # 4	4–Adult
LDA6M3	MacIntosh # 3	3–Small Adult
LDA6M2	MacIntosh # 2	2–Pediatric
LDA6M1	MacIntosh # 1	1–Neonate
LDA6HA	Hyperangulated "H"	Adult
LDA6C3	Channeled Regular	3–Adult
LDA6C2	Channeled Small	2–Small Adult
LDA6C1	Channeled Pediatric	1–Pediatric
LD A6C0	Channeled Infant	0–Neonate

- All 3 blade geometries
- From neonates to adults
- Made from ISCC Certified Recycled Polycarbonate
- One-piece handle and blade minimizes cross-contamination
- Low profile
- Hermetic

ML2025-482 05/2026 Rev. B