Airtraq +Plus

ALL-IN-ONE ECO VIDEO LARYNGOSCOPE

LEARN MORE

About

The Airtraq® +Plus Video Laryngoscope sets a new standard in intubation with its versatile and intuitive design. Compatible with three blade shapes and four blade sizes, it features a dual-axis swiveling screen for effortless use from any position. Its lightweight Wi-Fi camera enables video recording and live streaming across multiple devices. Thoughtfully designed with sustainability in mind, it incorporates recycled polycarbonate blades and ECO-friendly packaging, making it a smart, planet-conscious choice for advanced airway management.



Video Unit



WiFi

Advantages

- Easy to use for all skill levels.
- Fits three blade shapes: MacIntosh, Channeled, & Hyperangulated.
- Fits four blade sizes: Neonatal, Pediatric, Small Adult, and Adult.
- Recording capability to verify intubation.
- Single-use, ECO-friendly blades from recycled polycarbonate materials.
- Compact, portable, and affordable design.
- High-resolution optics with HD image quality and anti-fog coating.
- Dual-axis swiveling screen with screen flip functionality enables versatile intubation.
- Wi-Fi or wired video sharing for large displays.
- Blade attachment triggers camera activation (3 seconds).
- Water resistant, easy to clean, with quick recharge battery.
- Three Years Warranty for Camera & Max View Display


Fits 3 Blade Shapes and 4 Sizes



MacIntosh (4 Sizes)

Hyperangulated (1 Size)

Channeled (4 Sizes)



Neonatal Pediatric Small Adult Adult

Clinical Use

- High-resolution optics and anti-fog technology.
- Dual-axis swiveling screen for flexible intubation angles.
- IP67 - Water resistant, easy to clean, with quick-rechargable battery.
- ECO-friendly design with recycled polycarbonate blades and sustainable packaging.

*IP67 - Indicates a product is both dust-tight and protected against temporary immersion in water, typically up to 1 meter deep for 30 minutes.

5/2026 ML2025-472 Rev. D

RECYCLED FSC packaging

Exclusively Distributed in the USA by

Mercury Medical® Superior Solutions

Dedicated to delivering clinically differentiated critical care technology that saves lives throughout the world.

(800) 835-6633

www.mercurymed.com

Airtraq +Plus

ALL-IN-ONE ECO VIDEO LARYNGOSCOPE

LEARN MORE



Airtraq +Plus Kit	LDA600	LDA600B
Video Unit	✓	✓
Charging Dock	✓	✗



Charging Dock



Nylon Storage Case

Battery Charger with Cables	✓	✓
Video HDMI Cable	✓	✓
Data Cable (USB-C to USB-C)	✓	✓
Silicone Protective Cover	✓	✓
Nylon Storage Case	✓	✓





MacIntosh Blades

- Extra-stiff, thin-profile design with authentic geometry.
- Optimal light intensity for direct laryngoscopy.
- Consistent glottis image size for Mac #3 and #4.
- Durable one-piece blade with integrated handle.
- Hermetically sealed to prevent contamination.

MacIntosh Size #1 Neonate	10/Box	**LDA6M1**
MacIntosh Size #2 Pediatric	10/Box	**LDA6M2**
MacIntosh Size #3 Sm Adult	10/Box	**LDA6M3**
MacIntosh Size #4 Adult	10/Box	**LDA6M4**





Hyperangulated Blades

- Optimized hyper-curved geometry for precision.
- Thin, extra-stiff blade for stability.
- Durable one-piece design with integrated handle.
- Hermetically sealed to ensure no contamination.

Hyperangulated	10/Box	**LDA6HA**



Channeled Blades

- Thin profile with 90° geometry for reduced force.
- Integrated channel simplifies ETT insertion from any position.
- Durable one-piece design with built-in handle.
- Hermetically sealed to help prevent contamination.

Channeled Size #0 Neonate	10/Box	**LDA6C0**
Channeled Size #1 Pediatric	10/Box	**LDA6C1**
Channeled Size #2 Sm Adult	10/Box	**LDA6C2**
Channeled Size #3 Adult	10/Box	**LDA6C3**



OPTIONAL
Max View Display

Blades made from RECYCLED POLYCARBONATE

Designed By:
Prodol Meditec S.A.

Max View Display (Optional)

10.4" Max View LCD Display	**LDA300**

With Pole Mounting Clamp, VESA Screws, and Charging Dock Cradle with Mounting Bar. *Stand Not Included

Exclusively Distributed in the US by:
Mercury Medical®
Superior Solutions