# Airtraq

GUIDED VIDEO INTUBATION

LEARN MORE

## About

Airtraq™ simplifies video laryngoscopy with its ETT channel guide aiding in safety and in reducing intubation time. The optional lightweight Wi-Fi camera facilitates video recording and auto recording options.




## Advantages

- Simple to use and easy to learn.
- Makes intubation easy for non-expert and expert clinicians.
- Guided channel leads the endotracheal tube (ETT) through vocal cords.
- Minimizes hyperextension of the neck.
- Provides a stable grip while intubating.
- Single use blades.
- Compact, portable, and affordable.
- Fully integrated high-resolution optics.
- Blades automatic turn on when camera is mounted. (Adult Sizes Only)
- Anti-fog coating allows for a clear view of airway.
- Improved image quality.
- Allows an indirect visualization of the upper airway.
- Improves the success rate of tracheal intubation.

## AirTraq SP Color-Coded Blades




Airtraq™ Wi-Fi Camera Full Kit
#LDA390KIT



Infant (Gray)
2.5 - 3.5 ETT
#LDA041 (2/Pkg)



Pediatric (Purple)
4.0 - 5.5 ETT
#LDA031 (2/Pkg)



Adult / Small (Green)
6.0 - 7.5 ETT
#LDA021 (2/Pkg)



Adult / Regular (Blue)
7.0 - 8.5 ETT
#LDA011 (2/Pkg)



Nasal (Orange)
#LDA061
(2/Pkg)



Double Lumen (Yellow)
#LDA071
(2/Pkg)

## Clinical Use

- Disposable color-coded video laryngoscope handle blade system with ETT guide channel & eyepiece for direct view.
- 4 Blade sizes. (Infant 0, Pediatric 1, Adult 2, Adult 3)
- 2 Specialty blades: Nasotracheal Blade and Double Lumen.
- Airtraq's™ channel simplifies ETT insertion reducing intubation time.
- No stylet needed. Avoids potential injuries. Eliminates stylet costs.
- 90º shape minimizes hyperextension and reduces force required.
- Reduces hemodynamic reactivity.
- Facilitates intubation from any position. (Face to face, etc)
- Video recording capability.

Manufactured by: Prodol Meditec S.A.

Exclusively Distributed in the USA by

Mercury Medical®
Superior Solutions

*Dedicated to delivering clinically differentiated critical care technology that saves lives throughout the world.*

**(800) 835-6633**

www.mercurymed.com

07/2024 ML2024-464 Rev. B