FLOW-SAFE II®

DISPOSABLE CPAP SYSTEM

LEARN MORE

About

Flow-Safe II® is a disposable CPAP System with intergrated manometer for verifying pressure.

Advantages

- 15 LPM flow generates 10 cm CPAP pressure.
- Consistent CPAP pressure.
- Disposable manometer.
- Pressure pop-off valve. (25 cm H_2O)
- Deluxe mask with adjustable forehead pads.
- Nebulizer capable.
- Capnography capable.
- Filter capable.
- Fits in first-in bag.
- Mask clips for easy access to patient.
- Adapts to ET Tube.
- FiO_2 up to 100% depending on flow rate, the patient's tidal volume, and respiratory rate.

Flow-Safe II®
(5/Box)



Flow-Safe II®
With Child Mask
(5/Box)



Clinical Use

- Provides CPAP to spontaneously breathing patients.
- Pre-hospital or in hospital use.
- Just attach to standard flowmeters/Flow determines the CPAP pressure
- Easy set up & training.
- Ideal for wherever CPAP therapy is appropriate.
- Compact and easy to store (No extra parts.)
- Color-coded manometer verifies delivered CPAP pressure.
- Comfortable head harness with quick disconnect clips for easy access to the patient's face.
- Forehead support reduces pressure on the nose and easily adjusts.
- Yellow Flow/CPAP pressure chart on all devices.

Mercury Medical®
Superior Solutions

Over 55 Years of Medical Device Sales and Development

Dedicated to delivering clinically differentiated critical care technology that saves lives throughout the world.

(800) 835-6633

www.mercurymed.com

07/2024 ML2022-443 Rev. A