Mercury Medical

Mercury Medical's Latex-Free CPR Bags

Mercury Medical's CPR Bags are designed to accommodate the most urgent and unique situations of on-site EMS, Hospital I.C.U., E.R., O.R. and Recovery Room needs.



Adult CPR Bag

Patients Over 30 kgs

Maximum Rate	92 BPM		
Maximum Volume	900 - 1500 ml		
Pressure Relief Valve	(Optional)		
Oxygen Concentration			
LPM	3	5	10
Rate	12	12	12
Tidal Volume	500 ml	500 ml	500 ml
W/Reservoir	68%	95%	100%
W/O Reservoir	34%	47%	70%
Measured Volume			
Bag	1600 ml		
Reservoir	2600 ml		



Child CPR Bag

Patients Between 7 and 30 kgs

Maximum Rate	168 BPM
Maximum Volume	360 ml
Pressure Relief Valve	40 cm H_2O
Oxygen Concentration	
LPM	10
Rate	20
Tidal Volume	250 ml
W/Reservoir	100%
W/O Reservoir	85%
Measured Volume	
Bag	500 ml
Reservoir	2600 ml



Small Adult CPR Bag

Patients Over 30 kgs

Maximum Rate	110 BPM		
Maximum Volume	500 - 750 ml		
Pressure Relief Valve	(Optional)		
Oxygen Concentration			
LPM	3	5	10
Rate	12	12	12
Tidal Volume	500 ml	500 ml	500 ml
W/Reservoir	68 %	95 %	100 %
W/O Reservoir	34 %	47 %	70 %
Measured Volume			
Bag	1000 ml		
Reservoir	2600 ml		



Infant CPR Bag

Patients Under 7 kgs

Maximum Rate	180+ BPM
Maximum Volume	175 ml
Pressure Relief Valve	40 cm H_2O
Oxygen Concentration	
LPM	4
Rate	30
Tidal Volume	40 ml
W/Reservoir	100%
W/O Reservoir	85%
Measured Volume	
Bag	240 ml
Reservoir	600 ml

Latex-Free CPR Bags including Infant, Child, Small Adult and Adult

11300 - 49th Street North · Clearwater, Fllorida 33762-4807 · Telephone: 727.573.0088 · Toll Free: 800.237.6418

Mercury Medical™

Mercury Medical's Latex-Free CPR Bags

6/Box

Adult CPR Bags



ORDERING INFORMATION

Product Number	Cush. Mask	Cont. Mask	O_2 Resv.	Aero Resv.	Peep Valve	Pop-Off	Lock Clip	Omni-Link	Other CPR Parts	Manometer	Exhalation Filter
10-55000*		X	X								
10-55002*											
10-55003			X								
10-55004	X		X								
10-55007	X			X							
10-55008*				X							
10-55023*	X		X						1,2		
10-55024*	X			X	X						
10-55025*	X		X		X						
10-55027*	X		X					X			
10-55028*			X					X			
10-55045*	X		X			X	X		16		
10-55048*	X		X								X
10-55059*	X		X		X					X	
10-55061*	X			X				26			X
10-55068	X		X								
ADULT SIZE WITH EXPANDABLE AEROSOL RESERVOIR									*12/Box*		
10-55700									23,24		
10-55701	X								23,24		
10-55702					X				23,24		
10-55703*	X				X				23,24		
10-55704*	X					X	X		23,24		
10-55705*	X								23,24		X
10-55713*	X								23,24	X	

*Custom configurations require 10 business days to manufacture and a minimum order of 10 boxes per order.

Adult CPR Bags with StatCO$_2$® End Tidal CO$_2$ Detectors



Product Number	Cush. Mask	O_2 Resv.	Aero Resv.	Expandable Aerosol Reservoir	StatCO$_2$ End Tidal CO$_2$ Detector	6/Box	12/Box
10-55064		X			X	X	
10-55065	X	X			X	X	
10-55066			X		X	X	
10-55067	X		X		X	X	
10-55711				X	X		X
10-55712	X			X	X		X

11300 - 49th Street North · Clearwater, Fllorida 33762-4807 · Telephone: 727.573.0088 · Toll Free: 800.237.6418

Mercury Medical™

Mercury Medical's Latex-Free CPR Bags

6/Box

Small Adult CPR Bag



ORDERING INFORMATION

Product Number	Small Adult Cush. Mask	Cont. Mask	O_2 Resv.	Aero. Resv.	Peep Valve	Pop-Off	Lock Clip	Omni-Link	Other CPR Parts	Manometer
10-55901	X								23	
10-55907	X			X						
10-55904	X		X							
10-55909						X	X		2, 23	
10-55910			X						2	X

6/Box

Child CPR Bag



ORDERING INFORMATION

Product Number	Cush. Mask	Cont. Mask	O_2 Resv.	Aero. Resv.	Peep Valve	Pop-Off	Lock Clip	Omni-Link	Other CPR Parts	Manometer
10-55100*		X	X			X	X			
10-55102*						X	X			
10-55103			X			X	X			
10-55104	X		X			X	X			
10-55106				X		X	X			
10-55107	X			X		X	X			
10-55113*	X		X		X	X	X			
10-55115*				X		X	X		3	
10-55119*	X			X	X	X	X			
10-55122*	X		X			X	X		3	
10-55142*	X		X			X	X		3,5	
10-55150*	X		X		X	X	X		3	
10-55156*						X			3	
10-55163*	X		X			X			3	
10-55166*	X		X			X	X		1	
10-55167*	X		X			X				
10-55168*	X		X			X	X		3	X
10-55170*	X		X			X	X		3,5	X
10-55171*	X		X		X	X	X			X
10-55172*	X			X						
10-55173*	X		X			X	X		18	
10-55175*	X			X	X	X	X		3	X
10-55177*	X			X		X	X			X
10-55179*	X			X		X	X		3,5	X
10-55180*			X			X	X		3	X
10-55182*				X		X	X			
10-55183*	X		X			X	X		26	
10-55186*	X		X			X	X			X
10-55187*	X		X			X	X		26	
10-55188*				X		X	X		3	X
10-55189*	X			X		X	X		18	X
10-55190*			X			X	X			X

11300 - 49th Street North · Clearwater, Fllorida 33762-4807 · Telephone: 727.573.0088 · Toll Free: 800.237.6418

Mercury Medical™

Mercury Medical's Latex-Free CPR Bags

StatCO$_2$ End Tidal CO$_2$ Detector Foil Pack added to CPR drawstring bag.
The StatCO$_2$ has a 2-year expiration date from date of manufacture.

*Custom Configurations require 10 business days to manufacture and a minimum order of 10 boxes per order.

Child CPR Bag with StatCO$_2$® End Tidal CO$_2$ Detectors



ORDERING INFORMATION

Product Number	Cush. Mask	O_2 Resv.	Aero Resv.	StatCO$_2$ End Tidal CO$_2$ Detector	6/Box
10-55185	X	X		X	X



6/Box

Infant CPR Bag

Product Number	Cush. Mask	Cont. Mask	O_2 Resv.	Aero Resv.	Peep Valve	Pop-Off	Lock Clip	Omni-Link	Other CPR Parts	Manometer	Mini StatCO$_2$
~~10-55200*~~		~~X~~	~~X~~		~~X~~	~~X~~					
10-55203			X			X	X				
10-55204	X		X			X	X				
10-55217*	X		X			X	X		2,5		
10-55223				X		X	X				
10-55224				X		X	X		5,20		
10-55229	X			X		X	X				
10-55234	X		X			X	X		5		
10-55241*	X		X		X	X	X				
10-55249*	X			X		X	X		5		
10-55250*	X			X		X	X		5,13		
10-55253	X			X		X			5,13		
10-55257*	X		X			X	X		5,8,21		
10-55264	X		X			X			5		
10-55266*						X	X		Size, 01 mask		
10-55271*			X			X					
10-55278*	X			X		X	X		5,21		
10-55280*	X			X		X	X		5	X	
10-55282*	X			X		X	X		11	X	
10-55283*	X			X		X	X			X	
10-55284*	X		X			X	X			X	
10-55287*			X			X	X			X	
10-55288*	X		X		X	X	X		5		
10-55289*			X			X			25		
10-55290*				X		X	X		Size 00/01 Masks		
10-55291*				X		X	X			X	
10-55292*	X		X		X	X	X			X	
10-55293*	X			X		X				X	
10-55294*			X			X	X		Size 00/01 Masks	X	
10-55295*	X		X			X			5,25		
10-55296*	X		X			X	X		18	X	
10-55297*	X		X		X	X	X		5	X	
10-55298*	X		X			X				X	
10-55299*	X		X			X	X		5	X	
10-552100*	X			X	X	X	X		5	X	
10-552101*	X			X	X	X	X			X	
10-552102*	X		X			X	X		26		
10-552107*				X		X	X			X	
10-552112*	X		X			X	X		5	X	X
10-552113*				X	X	X	X		27		

11300 - 49th Street North · Clearwater, Fllorida 33762-4807 · Telephone: 727.573.0088 · Toll Free: 800.237.6418

Mercury Medical™

Mercury Medical's Latex-Free CPR Bags

OTHER CPR PARTS

1	Small Adult Mask
2	Child Mask
3	Infant Cushion Mask
5	Neonate Mask
6	Infant Patient Valve
7	Glued Oxygen Tubing
8	Glued Reservoir
9	Pediatric Mask
10	Standard Oxygen Tubing

11	Assembled except Aerosol Tubing
13	Removable Aerosol
15	Strapped O_2 Tube to Aerosol
16	Disconnect Wedge
18	Exhalation Filter
19	Bag Hanger
20	Hook Ring
21	Airway Pressure Line
23	Expandable Aerosol Reservoir
24	Packaged 12/case
25	Assembled except Oxygen Tubing
26	7' O_2 Tubing
27	Custom Zip Lock Bag with Handle

• All CPR Bag Patient Valves listed have built-in PEEP Adapters.
• All CPR Bag configurations come fully assembled.

CPR ACCESSORIES

10-55310	PEEP Adapter
10-55330	PEEP Valve
10-55315	PEEP Adapter & Valve Package
10-55320	Manometer Port
10-55340	Disconnect Wedge
10-55360	Omni-Link 12"
10-55361	Omni-Link II 6"
10-55350	Manometer (Single-Patient Use)
10-55351	Manometer with Tubing and Stand
10-55375	Exhalation Filter

KEY TO TERMS USED

O2 Resv.	**Oxygen Bag Reservoir**
Cush. Mask	**Cushion Mask**
Cont. Mask	**Contour Mask**
Aero. Resv.	**Aerosol Tubing Reservoir**

11300 - 49th Street North · Clearwater, Fllorida 33762-4807 · Telephone: 727.573.0088 · Toll Free: 800.237.6418