
Mercury Medical®

Your Need ... Our Innovation®

ONLY ONE



NIFometer
(to -40 cm H_2O)
10-340-00
20/Box



NIFometer
(to -60 cm H_2O)
10-340-02
20/Box

Single-Patient-Use

NIFometer®

NEGATIVE INSPIRATORY FORCE METER

The quick & convenient way to measure NIF.

QUICK AND EASY

- Small, compact, self-contained. Connects directly to ET tube eliminating need for separate tubing and connectors.
- Includes built-in occlusion button that you simply depress during test.

SINGLE PATIENT USE

- Helps eliminate risk of cross-contamination.

ACCURATE

- Accurate to ± 2 cm H_2O -20 to -30 cm H_2O (NIF 0-40) & ± 3 up to -15 cm H_2O & ± 5 from -15 to -60 cm H_2O (NIF 0-60).
- No calibration required.
- Includes memory indicator pointer that records highest force reached and can then be reset.

CONVENIENT

- Lightweight weighing only 30 Grams.
- Provides unobstructed view of patient during test.
- Small and compact, fits in pocket – readily available for added efficiency.
- Also available in kit with elbow, mouthpiece and nose clip.

11300 - 49th Street North • Clearwater, Florida 33762-4807 USA

Tel. (800) 835.6633 • (727) 573.0088 • Fax (800) 990.6375

www.mercurymed.com

Since 1963

over ►

Mercury Medical®

Your Need ... Our Innovation®

ONLY ONE



NIFometer
(to -40 cm H_2O)

NIFometer
(to -60 cm H_2O)



NIFometer Kit
Available with NIFometer
to -40 cm H_2O & to -60 cm H_2O

Single-Patient-Use

NIFometer®

NEGATIVE INSPIRATORY FORCE METER

The quick & convenient way to measure NIF.

SPECIFICATIONS

Range	**To -40 cm H_2O**
Variance	± 2 cm H_2O -20 to -30 cm H_2O
Fitting	15 mm I.D./22 mm O.D.
Weight	30 grams nominal

Range	**To -60 cm H_2O**
Variance	± 3 up to -15 cm H_2O & ± 5 from -15 to -60 cm H_2O
Fitting	15 mm I.D./22 mm O.D.
Weight	30 grams nominal

ORDERING INFORMATION

Part #	Description	Packaged
10-340-00	NIFometer (to -40 cm H_2O) Single-Patient-Use	20/Box
10-340-01	NIFometer Kit, (to -40 cm H_2O) Includes: elbow, mouthpiece & nose clip	10/Box
10-340-02	NIFometer (to -60 cm H_2O) Single-Patient-Use	20/Box
10-340-03	NIFometer Kit, (to -60 cm H_2O) Includes: elbow, mouthpiece & nose clip	10/Box

Patent # US 6,854,334 B2

11300 - 49th Street North • Clearwater, Florida 33762-4807 USA

Tel. (800) 835.6633 • (727) 573.0088 • Fax (800) 990.6375

www.mercurymed.com

© 7/2020

ML2012-220 Rev. B

Since 1963